



# Phase-1 Upgrade of the CMS Hadron Calorimeter Endcaps

Caleb J. Smith (Baylor University)
October 26, 2018



# CMS HCAL Phase 1 Upgrade






#### New photosensors

 Replace Hybrid Photo-Detectors (HPD) with Silicon Photomultipliers (SiPM)


#### New readout electronics

QIE cards with custom QIE ASIC chips






#### Phase 1 Readout Module (RM)



Silicon Photomultipliers (SiPMs) and Peltier





4.8 gbps optical data links

JTAG adapter board



#### Phase 1 Readout Box (RBX)









## HCAL Endcap (HE)

CMS


- 2 HCAL Endcaps (±z)
- 18 wedges per Endcap
  - 1 Readout Box (RBX) per wedge
  - 20 deg. (ф) per wedge









#### Motivation








integrated luminosity (fb<sup>-1</sup>)

Talk: HCAL HE Status by Pawel De Barbaro


- Radiation damage to Hybrid Photo-Detectors (HPD) occurred sooner and was more severe than expected.
- The HE HPDs were replaced with SiPMs in the Phase 1 upgrade in early 2018.








#### From HPDs to SiPMs





# HPD (Hybrid Photodiode) 150 GeV muons Sum of 12 layers 1.5 10 10 Energy [GeV]

#### **HPDs to SiPMs**


- x 400 higher gain
- x 2.5 PDE
- Lower voltage (8kV to 70V)
- Reduced noise



A.U. / GeV

Calorimeter response to muons from HE prototype wedge at H2 beamline







Phase-1 Upgrade of the CMS Hadron Calorimeter Endcaps



#### New Geometry





7.0 m

11.15 m -


- Increase in number of channels.
- Increase in number of depths (each color is a depth).
- Provides better longitudinal segmentation.




#### Installation

#### Procedure

- Remove existing electronics.
- Install new electronics.
- Test control and data links.
- Take FiberID, LED, and Laser data.
- Perform radiation sourcing scan with Co60.
- Calibrate using Co60 data.











#### Source Driver

- Radioactive source (Co60) is pushed with a wire through tubes in the detector.
- As the wire extends, the source passes near different scintillator tiles.
- Response measured in every channel.
- Used to verify end-toend channel mapping and cable connections.
- Used for startup calibration of detector.











Wire extension (mm)









Improved uniformity vs. phi after the upgrade!








(Charge - Mean) / Mean

- Comparison of 2017 sourcing (HPD) and 2018 sourcing (SiPM).
- Improved uniformity in 2018.




#### Conclusion



- Installation and commissioning were completed on schedule during the 2017–18 year-end technical stop.
- The HE Phase 1 upgrade brings improved detector response and uniformity.
- The HCAL Barrel (HB) Phase 1 upgrade will be installed during Long Shutdown 2 (2019–2020).


Thank you for your time!






# Backup









- http://feynman.physics.uiowa.edu/HESourcing18/Others/SourceResponse HEP.pdf
- http://feynman.physics.uiowa.edu/HESourcing18/Others/SourceResponse HEP 2017.pdf
- http://feynman.physics.uiowa.edu/HESourcing18/Others/SourceResponse HEM.pdf
- <a href="http://feynman.physics.uiowa.edu/HESourcing18/Others/SourceResponse">http://feynman.physics.uiowa.edu/HESourcing18/Others/SourceResponse</a> HEM 2017.pdf








Run308040\_HEP15\_PHI55\_LAYER13\_SRCTUBE\_S\_eta17\_phi55\_depth3



 Megatile cable is connected backwards




Reel position (mm)


**O ?** 




# Commissioning with LED




- Low intensity LED runs were used to tune the SiPM gains.
- The bias voltages are selected such that the gains are tuned to 40 fC (charge measured for one electron).








## Commissioning with Laser



- Laser to Megatile
  - Only measured by a subset of channels.

- Laser to Calibration Unit (CU)
  - HEP receives more light than HEM due to laser splitting.









#### HCAL endcap sectors 15/16

Talk: CMS Status Report by Jim Hirschauer

- Following power interlock on June 30, two endcap sectors are not functional.
  - 40° in one endcap, 2% of HCAL coverage.

- Five-week campaign led to full understanding:
  - On power up after interlock, 10V power supply
     (PS) unable to read internal calibration.
  - PS sent 22V/10ms pulse to detector
    - exceeded its own 14V max rating
    - damaged on-detector components with 12V rating.












- HCAL installed secondary safety system to mitigate risk of damage from potential future transients.
- PS manufacturer working to understand and address
  - why the PS fails to read its calibration
  - why the PS sends high voltage when the calibration fails
- Physics impact:
  - trigger rates are OK
  - effect on MET resolution is small but measureable
  - PF reconstruction reduces impact of loss.
  - Additional modifications of reconstruction in progress to minimize impact.


Talk: CMS Status Report by Jim Hirschauer



#### HCAL Endcap (HE)





- HCAL Endcap (HE)
- 18 wedges
  - 1 RBX per wedge
  - 20 deg. per wedge
- Colors correspond to power supplies
  - 1 or 2 RBX per power supply






#### HCAL DQM









#### Electrons (n)












#### Electrons (ф, рт)












#### Jets (n)












#### Jets (ф, рт)









#### MET and H<sub>T</sub>

