



- drifting charges are extracted from the liquid to the vapour phase by an electric field in the liquid of around 2 kV/cm
- the charges once in pure argon vapour are multiplied inside 50x50 cm2 area LEMs (Large Electron Multipliers) assembled side by side
- the amplified charges are collected on a 2D segmented anode




3x1x1 construction 2015-2016













Mar 2017 - Operation 'frozen' due to cryostat issues






June 12th Recirculation started June 15th - evidence of extraction from LAr to GAr


June 21st 2017 - First tracks observed

Data taking July-November 2017





hadronic and EM showers from cosmic rays recorder in the 3x1x1



~400 k cosmic events recorded.

Article published in 2018



Journal of Instrumentation

Journal of Instrumentation, Volume 13, November 2018

OPEN ACCESS

A 4 tonne demonstrator for large-scale dual-phase liquid argon time projection chambers

B. Aimardⁱ, Ch. Alt^a, J. Asaadi^t, M. Auger^b, V. Aushev^v, D. Autiero^c, M.M. Badoi^r, A. Balaceanu^r, G. Balikⁱ, L. Balleyguier^c +Show full author list
Published 6 November 2018 • © 2018 CERN

- Submitted to Arxiv in June 2018,
 published in JINST November 2018
- <u>Detailed summary of the detector design</u> and performance in 60 page article.

1	Co	Contents						
2	1	Introduction	3					
3 4 5 6	2	Overview of the experimental apparatus 2.1 Experimental setup	4 6 10					
7 8 9 10	3	Cryostat and cryogenic system 3.1 The cryostat	11 11 14 14 16					
12 13 14 15 16 17	4	Description of the detector 4.1 TPC drift cage and drift field high voltage	17 17 19 19 22 24 26					
19 20 21 22 23 24	5	Ancillary instrumentation and detector control system 5.1 Pressure and temperature monitoring	28 30 30 32 33 34					
25 26 27 28 29 30 31 32	6	Detector readout, data acquisition and processing 6.1 Charge readout overview	36 37 39 40 43 44					
33 34 35 36 37 38	8	Detector performance 7.1 Cryogenic system performance and stability of the liquid argon surface 7.2 Observation of prompt scintillation and electroluminescence 7.3 TPC charge readout and operational feedback	47 47 48 51 51 54					

3x1x1 - key technological milestones





- Performance of the cryogenic system and of the 23m³ membrane cryostat.
- Extraction of the ionisation charge over an area of 3m².
- Amplification in pure argon vapour by a combined operation of multiple 50 × 50 cm² LEMs
- Readout of the signal on two collection planes with strips of up to 3 m length and performance of analogue front-end electronics.
- Detection of prompt and secondary scintillation.

Cryostat and cryogenics



The goal of the cryostat and cryogenic system is to ensure during extended operating periods:

1. stable thermodynamic conditions

- Uniform and stable gas argon density near 1 Atm and 90 K completely isolated from outside
 P,T fluctuations
- operation of Front End electronics near 110 K
- compensate the heat input of cryostat by active cooling

2. stable liquid level during liquid recirculation


- Avoid slow fluctuations of level which might be induced by regulation of cryo system
- avoid fast fluctuations (=waves) which might be induce from local heat inputs (i.e bubbling) and defect in cryostat insulation.

3. high liquid argon purity:

• better than 3 ms drift or less than 100 ppt 02 equivalent impurities


Remark on the liquid argon level






an important point on requirement of level position:


- for a given $\Delta V_{LEM-grid}$ the extraction *field* depends on the position of the LAr level.
- At sufficiently large $\Delta V_{LEM-grid}$ (>~2.5 kV) the extraction efficiency (hence the gain) is near maximal and therefore almost independent of the liquid level.
- The boundary conditions are (of course) that the liquid should not touch the LEMs and the grid stays immersed.
- ->LAr level and CRP alignment to the LAr level within 1 mm and monitoring within 100 microns






Cryostat and cryogenics - liquid argon level











Cryostat and cryogenics - performance



- Temperature is uniform across the CRP.
- measured gradient near LAr surface is uniform at ~2 K/cm


 Pressure is stable at 999±1mbar decoupled from external pressure variations







Charge amplification and readout - the CRP









Precisely machined frame that allows to position the extraction-grid, LEMs and anodes with well determined inter-stage distances.

The frame is remotely adjustable to the liquid argon level with a precision of around 100 microns

<u>Demonstrate</u>:


- ⇒collective operation with gain of twelve LEMs arranged side by side
- ⇒extraction on 3m2 area and charge readout on 3m long strips

Charge amplification and readout - the CRP




LEMs and anode

- <u>large area LEMs (50x50 cm2)</u> o(150) holes per cm^{2.} Intrinsic geometrical parameters (rim, hole, pitch) optimised from test on smaller devices.
- Two view collection anodes with fine 3 mm segmentation: designed to provide identical charge sharing between both views with strip capacitance pe unit length as low as 150 pF/m



literature NIM A617 (2010) p188-192 NIM A641 (2011) p 48-57 JINST 7 (2012) P08026 JINST 8 (2013) P04012 JINST 9 (2014) P03017 JINST 10 (2015) P03017






• CRP Mechanical behaviour at cold: open bath LN2 test (shock cooling) with mechanical deformations monitored by photogrammetry



Front end amplifiers







DP detector design that allows:

- the operation of FE amplifiers at cryogenic temperatures (stable ~140 K as measured during 3x1x1 data taking)
- the positioning of the FE amplifiers close to the anodes.
- the **accessibility** of the FE cards at any time during detector operation for replacement or maintenance.

Double slope gain with dynamic range up to 1200 fC



Equivalent Noise charge vs detector capacitance



main argon volume


FE ASIC amplifiers located

feedthroughs isolated from

inside custom made "signal"

Electronic noise






- Noise at cryogenic temperature stable at around 1550 e⁻¹
- M.I.P Signal on one channel expected around (infinite electron lifetime):
 - ~10 ke for effective gain 2 (S/N=10)
 - ~100 ke for effective gain 20 (S/N=100)

✓ Noise is below 2'000 electrons on 3 m long strips

Cathode VHV during data taking





stable:


- current of 50 ±0.08 uA (as expected from the divider)
- voltage 56.5 ± 0.01 kV

max voltages:

- -on feedthrough alone: -300 kV JINST 12 P03021 arXiv:1611.02085
- -on 3x1x1 cathode -65 kV for the entire duration of the data taking (no discharges, stable current)

CRP HV during data taking





- The extraction grid discharged at around -5 kV while the nominal voltage (for an effective gain of 20) is set to around -6.8 kV.
- Evidence (after dedicated tests and visual inspection) points to a group of 32 wires which were not correctly tensioned and a faulty HV contact.
- As a consequence we were not able to operate the detector at its full potential with the nominal voltages shown above.
- Nevertheless we collected a large amount of good quality data at various field configurations.

CRP HV scanned ranges and data taking settings



Table 9. Electric field settings explored during and data taking with the $3 \times 1 \times 1$ m³ TPC (see text for more details). The values quoted as nominal correspond to G_{eff}^{∞} of 20 as achieved in the dual-phase TPC with the 10×10 cm² readout area [15, 22].

Electric field (kV/cm)	Scanned range	Data taking settings	Nominal
Induction	1–5	1.5	5
Amplification ³³	23–31	28	33
Extraction (in liquid)	0.6–2.5	1.7	2
Drift	0.18-0.7	0.5	0.5

scanned ranges with max electric fields obtained independently on each element by reducing the other field values in order to be compatible with the overall limitation imposed by the grid HV.


Stable data taking settings. best compromise taking into account the limitation on the absolute grid HV

Data taking July-November 2017




more than 400 k cosmic events collected

HV: data taking settings



Crossing muons





Crossing muons (MIPs): easy to reconstruct in 3D and deposits a know amount of charge of ~10 fC/cm in liquid argon (at 500 V/cm and w/o attenuation due to drift) =>characterise the TPC in terms of amplification (effective gain), LAr purity, charge sharing between views,...

The signal measured on one view

charge sharing between views **50% charge sharing**

effective gain goal 20

10 fC/cm (at 500 V/cm drift)

$$<\frac{dQ}{ds}>_{view_i} = F_{share}(view_i) \times \varepsilon_{extraction} \times G_{LEM} \times \varepsilon_{induction} \times <\frac{dQ}{ds}>_{MIP}$$


Charge sharing between views



anode requirements:

- low capacitance per unit length (<2000 enoise equivalent)
- 50:50 charge sharing between the strips. Tested at the cm strip scale (JINST 9 (2014) P03017)

50:50 charge sharing demonstrated here on 3x1 m2 strips





JINST 9 (2014) P03017





From 3D reconstructed tracks: calculate the $\Delta Q/\Delta x$ per view






charge is equally shared between views

Liquid Argon Purity







Select tracks that cross the entire 1 m drift.
Compute the dQ/dx as a function of drift.

$$dQ/ds \propto e^{-t_{drift}/\tau_e}$$


✓electron lifetime of few milliseconds achieved





Charge collection - waveforms






√Two views with similar unipolar response, signals well above the ambient noise.

 similar amplitudes on both views: visual confirmation 50/50 charge sharing between views

Effective gain




increase of the mean and MP value of the dQ/dx distributions as a function of the LEM fields. 28 kV/cm of LEM field corresponds to an operation of the 3x1x1 TPC at $G_{eff} \approx 3$.





Scintillation in liquid and gas argon





- Clearly visible light from primary and secondary scintillation.
- First check of correct event matching between charge and light events. Correlation between the quantity of light and charge detected between matched events.
- Comparison between the light simulation and data.





Conclusions



- Despite technical issues that prevented from stably operating the TPC at G_{eff} of 20, we have demonstrated the high quality of the dual-phase TPC imaging and provided a proof of principle that such devices can be constructed and operated at the ton scale.
- A total of about 400 k high quality cosmic events retrieved in many different field configuration.
- Good performance of cryostat and cryogenic system: purity compatible with ms electron lifetime and LAr level stable at sub-mm scale.
- Stable 500 V/cm drift field over one meter
- First time ever, charge extraction over 3m² area and LEM amplification demonstrated on 50x50 cm² modules.
- stable performance of FE electronics, low noise and no evidence of channels dying during detector operation.

Conclusions II




- First LAr TPC operation in a membrane tank and excellent performance of cryogenic system
- proof of stable charge extraction over a 1 cm gap on a 3 m² surface in 17m3 of LAr
- First time operation of multiple side by side 50x50 cm² LEMs with gain
- Operation of an entire 3m2 CRP (anode, LEM + grid) and gained a detailed understanding of the capacitive couplings between each stage.
- Many methods developed to monitor level (precision level meters, LEM-grid capacitive measurements ..)
- First time use in a LAr TPC of accessible cold front end electronics: they have shown to be robust to discharges and offer excellent noise performance even with readouts of ~500 pF (3 m strips)
- Stable operation of the drift cage with cathode at 60 kV
- excellent quality of charge readout offered by two collection views and gain.
- Noise filtering and track reconstruction, many complex "neutrino like" events, unique opportunity to test and tune reconstruction algorithms.
- light data: observation of S2 at different extraction and LEM fields, effect of purity, drift field,... on S1. In general comparison of response with simulations to validate and tune photon propagation models for future LAr-TPCs



Gas Argon evolution purity during purge




Evolution of the impurities measured in the gas during open, closed loop piston purge and cooling down. The measured temperatures inside the gas near the bottom, middle and top of the main cryostat volume are indicated in blue.




Front end electronics





temperature FE: ~140 ±5 K without SGFT cooling coils. (see Franco's slides for effect of cooling)




4000

2000

Charge amplification and readout - the CRP





The system can be represented by a network of capacitively coupled electrodes including the anode, the LEM down, the LEM up, the grid and the field cage. We have measured all the involved capacitance and define a ramping up procedure where all the currents involved are adjusted taking into account the capacitive couplings.

Table 4. Approximate inter-electrode capacitance measured in the $3 \times 1 \times 1$ m³ (see text for nomenclature).

Electrode	distance (mm)	capacitance
		in GAr (pF)
3 m (1 m) anode strip to GND	-	480 (160)
$50 \times 50 \mathrm{cm}^2$ LEM-top to anode plane	2	1000
$50 \times 50 \mathrm{cm^2}$ LEM-top to LEM-bot	1	7000
$50 \times 50 \mathrm{cm^2}$ LEM-top to neighbour LEM-top	4.5	≤1
$50 \times 50 \mathrm{cm}^2$ LEM-bot to extraction grid	10	150
$3 \times 1 \text{ m}^2$ cathode to GND	200	4000






Figure 30. Left: comparison of noise measurements at warm (red points) and cold (black points) with the slow control cables disconnected. Right: comparison of noise measurements taken at cold by progressively reconnecting various equipment.







Figure 19. Reading from one CRP level meter (left) and $50 \times 50 \,\mathrm{cm^2}$ LEM to extraction grid capacitance (right) as a function of the vertical position of the CRP given by the motorised system.


Dual phase R&D prior to the 3x1x1 2005-2015





 $\Delta Q_0^{200} / \Delta s_0$ [fC/cm]


CRP: optimal design of LEMs, anodes, extraction grid, FE amplifiers
















all components designed to be scalable to k-tonne devices