




LArIAT in 12+3 Minutes



What is LArIAT?

- Liquid Argon In A Testbeam
- LArIAT is a Liquid Argon Time Projection Chamber (LArTPC) in a beam of **known** particles.


Cryostat and
TPC
refurbished
from the
ArgoNeuT
experiment



[1]


LArTPC Technology

- 1) Charged particles ionize and excite argon.
- 2) Argon scintillates, and the PMTs record the photon within nanoseconds, marking initial time for the event.
- 3) Ionization electrons drift towards induction and collection planes on the order of hundreds of microseconds.
- 4) Electrons induce current in the planes and land on the collection plane.




[1]

Pions and Protons!





Track Reconstruction




Why do we care?

- LArTPCs are the neutrino detectors of choice for the SBN program and DUNE
- Using calorimetry and topology to distinguish E&M showers between **photons** and **electrons**
- LArIAT **calibrates** future neutrino experiments
 - Measurements of how neutrino interaction products behave in liquid argon
 - Methods to explore liquid argon technology




Test Beam Facility

- LArIAT ran in Fermilab's Test Beam Facility from 2015 to 2017
- Tertiary beam; Particles with momentum between 0.3 GeV/c and 1.4 GeV/c



[2]


Beamline



- 13° between secondary and tertiary beam.
- **Wire Chambers** record timing and location info.
- **Magnets**: Run at 60A and 100A. Curve particle path.
- **TOF** system made of PMT scintillator paddles.
 - Upstream: 6 cm x 10 cm; Downstream: 14 cm x 14 cm


LArIAT Analysis

- Wire Chamber Momentum Reconstruction
- Wire Chamber Particle Identification (PID)



LArIAT Analysis

- Wire Chamber Momentum Reconstruction
- Wire Chamber PID
- Wire Chamber to TPC matching




π^- and K^+ Cross-Section

Why?

- Novel Measurement
- Affects detection and measuring energy of hadrons in a neutrino experiment


(π^-, Ar) Total Hadronic XS

LArIAT Preliminary



(K^+, Ar) Total Hadronic XS

LArIAT Preliminary




[3]

Pion Absorption Study

- Exploration of absorption π - kinematics
 - # protons
 - Energy of emitted protons
 - Energy loss

Why?

- Understanding the final state of neutrino interactions, where we don't necessarily see the pions.
- What can pions in the nucleus do?





Antiproton Annihilation

- Study the hadron star topology of proton-antiproton annihilation at rest.

Why?

- Calibrate DUNE for neutron-antineutron oscillations, which also occur at rest.




Michel Electron Calorimetry

- Novel method of analysis!
- Using liquid argon scintillation light to reconstruct energy of 5-50 MeV electrons from cosmic muons.

Why?

- Better understanding of low energy electron showers, such as from solar or supernova neutrinos in DUNE.



Summary

- 1) LArIAT is made up of the **beamline** and **LArTPC**
- 2) **Test beam** contains particles present in a neutrino interaction (muons, pions, electrons, kaons, protons)
- 3) **Beamline** has PID capabilities
- 4) **LArTPC** provides calorimetric and topological information
- 5) We use this to understand the behavior of these particles in liquid argon to **calibrate** LArTPCs for experiments like MicroBooNE, SBND, DUNE.
- 6) LArIAT allows us to develop new methods of analysis and make novel measurements!

References

- [1] LArIAT Detector Paper (JINST): LArIAT DocDB 4000
- [2] <https://neutrino.physics.fnal.gov/facilities/charged-particle-beam-facilities/>
- [3] NuInt 18 - 12th International Workshop on Neutrino-Nucleus Interactions in the Few-GeV Region Contribution
- [4] Michel Electron Calorimetry Paper: LArIAT DocDB 3646
- [5] Matt King Updated Beamline Reco: LArIAT DocDB 3965
- [6] <https://www.fnal.gov/pub/science/particle-accelerators/images/accel-complex-animation.gif>
- [7] Vincent Basque Antiproton Annihilation Presentation: LArIAT DocDB 4350

Backup Slides

LArIAT Analysis

- Monte Carlo Truth to develop measurement techniques

MC Truth




Data




- Box mimics physics but doesn't necessarily have realistic values
 - Create method for determining what's in the box (measurement).
 - Check if your method (which doesn't assume knowledge of contents of box) matches what's in the box.
- Apply method to data to make measurement.

LArIAT Analysis

- MC Truth to develop reconstruction algorithm
- Wire Chamber Momentum Reconstruction




The Test Beam



[6]

LArTPC Technology




- TPC Active Volume: (47 cm x 40 cm x 90 cm)
- Induction Plane and Collection Plane: 40 cm x 90 cm in Run I, II, IIIA. 30 cm x 70 cm Run IIIB.
- 3mm wire spacing.


Liquid argon temperatures and electric field measurements via the ACP track method

Run Period	Temperature (K)	Electric field (V/cm)
Run-I	89.4 ± 0.4	463.7 ± 10.6
Run-II	90.3 ± 0.3	465.9 ± 7.8
Run-III	90.3 ± 0.3	485.7 ± 12.0

[1]



Beamline



- 13° between secondary and tertiary beam. 10° bend.
- **Wire Chambers** record timing and location info.
- **Magnets**: gap height of 14.224 cm, gap width of 31.75 cm. Run at 60A and 100A.
- **TOF** system made of PMT scintillator paddles
 - Upstream: 6 cm x 10 cm; Downstream: 14 cm x 14 cm