



# Glass Thick-GEM

Takeshi Fujiwara (AIST), Yuki Mitsuya (Univ. Tokyo), Takashi Fushie (Radiment Lab. Inc.)

## Outline

1. Glass GEM fabrication process
2. Comparison with PEG3 and PEG3C
3. Applications
  - a. X-ray imaging with optical readout
  - b. Hadron Therapy
  - c. Neutron Bragg-Edge Imaging
4. Summary



1. Introduction
2. Glass gas electron multiplier with PEG3/PEG3C
3. Basic characteristics
4. X-ray imaging with analogue readout system

**Motivation: To develop a portable & versatile imaging detector.**

→ Sealed operation is essential

**Why sealed operation ?**

**1. Use of rare and expensive gases**

→  $^3\text{He}$  for neutron detection, Xenon for higher energy X-ray.

**2. High pressure operation ( or Low )**

→ Higher detection rate, smaller diffusion of electrons.

**3. Make a compact detector system**



→ Remove gas flow components, gas cylinder, purification system

**Non-outgassing MPGD is ideal**

→ Glass gas electron multiplier

# 1. Background & Motivation - X-ray imaging with Glass GEM


4




Cross section of the Glass GEM's hole

- ▶ Thick GEM made with glass substrate<sup>[1, 2]</sup>
- ▶ Why Glass GEM?
  - ▷ **Robust** – tolerant against discharges
  - ▷ **Rigid** – self-supporting structure, easy to handle
  - ▷ **Wet-etching** – uniform & Cylindrical Hole
  - ▷ **High gas gain** – up to 90,000 with single Glass GEM<sup>[3]</sup>
  - ▷ **High spatial resolution** – minimize charge spread


## An issue...



Metals inside via causes sparks



Micrograph of Glass GEM in early days



Discharges ruins the electrodes

[1] T. Fujiwara, et al., MPGD2011

[2] H. Takahashi, et al., NIM A, vol. 724, pp. 1–4, (2013)

[3] T. Fujiwara, et al., JINST, vol. 9, pp. 11007 - 11007, (2014)



## 2. Glass GEM fabrication process at AIST

5



UV exposure machine



Etching machine



DC sputtering machine for Cr




DC sputtering machine for Cu



Polishing machine




Plating machine



Laser microscope

Approximately 80% of the fabrication process are done in AIST (Many thanks to T. Fushie)

## 2. Glass GEM fabrication process




Optimizing the time and temperature

### 9. Ultrasonic bath


## 2. The new Glass GEM fabrication process at AIST

7

Old process




The New process



190μm

**Smooth electrode:** Uniformity of the electric field improves, and the GEM's stability improves








## 2. Comparison (PEG3 and PEG3C G-GEM)

Parameter comparison - Kapton, PEG3, PEG3C

	Unit	Kapton	PEG3	PEG3C
Volume Resistivity	$\Omega \cdot \text{cm}$	$\sim 10^{18}$	$8.5 \times 10^{12}$	$4.5 \times 10^{14}$
Bending Stress	MPa		>65	>150
Young's modulus	Gpa	18.6	79.7	90.3
Relative Permittivity		3.55	6.28	5.26
Heat conductivity	W/mK	$\sim 0.3$	0.795	2.72
Thickness	$\mu\text{m}$		400 – 800	100 - 700

Glass type



PEG3 G-GEM


### PEG3

**Low resistivity ( $8.5 \times 10^{12} \Omega\text{cm}$ )**

**= Less charge-up on the glass substrate**

= Stable operation is expected even with high rate.

Ceramic type



PEG3C G-GEM


### PEG3C=Crystallized PEG3

- Higher resistivity than PEG3 (may lead to more charge-up.)
- **Higher mechanical robustness increases available detector designs.**  
**e.g. thinner substrate** ( $\sim 100 \mu\text{m}$ )  $\rightarrow$  Thinner substrate = Less charge-up.
- Chemically stable, more phobic to moisture than PEG3.


## 2. G-GEM with PEG3/PEG3C

10

PEG3





PEG3C



Glass type


Ceramic type




**G-GEM Photos ( 10×10 cm<sup>2</sup> )**

## 2. G-GEM with PEG3/PEG3C


11




Max gain  $3.3 \times 10^4$



22.2% FWHM @Gain~10<sup>4</sup>



Max gain  $2.5 \times 10^4$




25.7% FWHM @Gain~10<sup>4</sup>

**PEG3 has higher gas gain and better energy resolution**


## 2. G-GEM with PEG3/PEG3C

12



**PEG3**

PEG3 was stable during high count rate measurement



**PEG3C**



Gain decrease due to the charge-up




**X-ray imaging with optical readout**

### 3. Optical readout Glass GEM<sup>[8,9]</sup>

14



- ▶ Glass GEM × scintillation gas × mirror × optical camera
- ▶ Convert radiation into visible light
- ▶ Optical mirror to prevent the CMOS from irradiated directly with X-rays
- ▶ ITO coated transparent electrode were used as an anode & optical window



[8] T. Fujiwara, et al., JINST, Vol. 8, No. 7 (2013)

[9] T. Fujiwara, et al., NIM A, 850 (2017)

### 3. Optical readout Glass GEM (X-ray imaging setup)

15

*DAQ and control PC*

*CMOS Camera*

10 sec integration

*dark box*

*Scintillating gas filled chamber*

*Glass GEM  
(100×100mm)*

*X-ray tube  
20 kV, 100  $\mu$ A*

*Object (a hornet)*

*X-ray*

*mirror*

*Rotating stage*

*Gas scintillation*








**CMOS Camera**  
Hamamatsu  
ORCA-Flash4.0 V3



**X-ray tube**  
Hamamatsu  
L9631

### 3. Optical readout Glass GEM - X-ray imaging


16





### 3. Optical readout Glass GEM - 3D X-ray Computed Tomography (CT)


17



**Scanned object (a hornet)**



**360 rotated X-ray images (1 degree step)**



**Reconstructed 3D CT images**

High resolution and high contrast X-ray 3D CT images of a hornet which consists of soft tissues were successfully reconstructed.



# **Dose Imager for Hadron Therapy**

# 1. Background and motivation – what is the issue for dosimetry?

19

## Conventional detector 1: Ion chamber

- ▶ Standard in clinical use
- ▶ Sharp Bragg peak is achievable
- ▶ Peak-to-Plateau ratio up to 4~5.
- ▶ **Spatial resolution is not enough (5mm)**
- ▶ **Takes time for each measurement**



Ion Chamber array  
(5mm pixel)


Peak-to-Plateau ratio

4.2 ~ 5.0

Proton/Carbon beam

## Conventional detector 2: Solid detectors

- ▶ Great spatial resolution (films, imaging plates, scintillator screens, semi-conductors)
- ▶ Radiation hardness would be an issue
- ▶ **Saturated in Bragg Peak**



Peak-to-Plateau ratio

1.3 ~ 3.2

Proton/Carbon beam


**Solid detectors has a quenching effect in high-LET radiation.**

- Energy deposition density at Bragg Peak is larger than the density of luminescence center.

**MPGDs have very little quenching effect and high spatial resolution**


# Detector construction (side view)

20







Updated to Brighter lens, high sensitive camera and shorter camera mount.






The pixel value of obtained image from CCD camera is plotted for each depth.






- Bragg curve of 290 MeV/u Carbon obtained with pixel values of CCD camera coupled with a scintillating gas-filled G-GEM chamber.
- **Peak-to-Plateau ratio was 4.48 (Ion-chamber), 4.51 (Glass GEM). (100.1%)**

# **Bragg-Edge Neutron Imaging Detector**




Reaction possibility -> Low  
Self absorption -> Low




Reaction possibility -> High  
Self absorption -> High


- ▶ B-10 are low price neutron converter
- ▶ Detect the charge from ionized  $\alpha/\text{Li}$  particle
- ▶ Self absorption of the particle is an issue
- ▶ Charged particle cannot escape to counting gas








- ▶ Slanting the absorber layer towards the incoming beam
- ▶ Grazing incident angle allows a larger proportion of neutrons to be absorbed in the first few microns of the layer
- ▶ It results secondary particles have a higher probability of escaping into the counting gas.
- ▶ This leads to increase neutron detection efficiency.







Micrograph of converter





3-D image






Picture of the micro-structure converter



Glass GEM








Spatial resolution = 0.8mm (FWHM)

# First Bragg-edge spectrum obtained with Glass GEM in AIST (Preliminary)

30






- ▶ Introduced the new fabrication process of Glass GEM
- ▶ **Application of Glass GEM detector**
  - ▷ X-ray imaging
  - ▷ Dose imaging detector for hadron therapy
  - ▷ Neutron Bragg-edge imaging detector
- ▶ Glass GEM is now open to everyone, and collaborators are always welcome

**Thank you for your kind attention.**

# Backup Slides


### 3. Initial experiment : $^{55}\text{Fe}$ (5.9keV X-rays) & PMT [7]

33




Scintillation light from Glass GEM is proportional to charge

Charge and light measured at same time



Good match




Emits 85,000 phtons/keV @gas gain 9,000


# 1. Background and motivation – what is done in hadron therapy?

34

## Hadron therapy



## Quality assurance



Before the treatment, **quality assurance of treatment** is done with **precise dose measurement**.

Personalized dose is measured and must be confirmed that has good enough agreement with the treatment plan.