#### Status and Performance of the Multi-anode MCP-PMT for the Belle II TOP counter

K. Matsuoka (KMI, Nagoya Univ.)
S. Hirose, T. Yonekura,
T. Iijima, K. Inami, R. Mizuno, Y. Kato, Y. Sato,
D. Furumura, K. Suzuki, T. Hayakawa

Micro-Channel Plate Based Detectors, Argonne, Dec. 2, 2014



Mass production of the MCP-PMT for Belle II TOP

- Performance of the MCP-PMTs
- Lifetime of the MCP-PMT photocathode
- Summary


# Mass production of the MCP-PMT

- The MCP-PMT for the TOP counter is a newly developed photon sensor in collaboration with HAMAMATSU.
   First time for the mass production of 512+α MCP-PMTs.
  - Production by HAMAMATSU and performance check at Nagoya are like the two wheels of a cart.
    - Feedback of our measurement is essential to successful production of MCP-PMT of good quality.






## **MCP-PMT** specification




Use 32 PMTs per TOP module, 512 PMTs in total.

# **Requirements for the MCP-PMT**



# Strategy of the mass production

6



These do not depend on HV.

## **Progress of the mass production**


7

Mass production started in March 2011 and finished in March 2014.



#### **Performance of the MCP-PMTs**


### QE test setup




#### **QE** measurement

10

Scan the photocathode at 18 x 18 points x 20 wavelengths.




# QE (bad case)



The QE decrease at the corners indicated a defect of the vacuum sealing between the window and the tube.  $\rightarrow$  Shored up the sealing.

12 **QE** 




Average QE at the peak  $\lambda$  > 28%

#### Laser test setup

13


Irradiate single photons to each channel one by one. Take data at 7 different HVs.

Automated (Just leave it for 2 h after click start button).




### Laser test readout electronics


Special variable-gain amp to keep ~100% readout efficiency even at a very low MCP gain below  $5 \times 10^5$ .




# **Calibration of the readout electronics**




#### **Gain measurement**



# <sup>17</sup> Gain uniformity



### TTS measurement




## **Collection efficiency measurement**

 $\rm CE \propto Number \ of \ TDC \ hits \ / \ QE$ 

19

Variation of the laser intensity is corrected by the reference PMT.




Higher CE of the coated MCP by ~17% than the conventional one.


#### Laser test setup in 1.5 T

20

Fix the MCP-PMT not to be moved by magnetic force.Use an MPPC (insensitive to magnetic filed) as a reference.



Take two data sets with the magnet off and on.




TTS(1.5 T)

CE(1.5 T)


Gain(1.5 T)

#### The gain decreases in 1.5 T.




Rate of decrease is different between the conventional and coated MCP-PMTs.

#### <sup>22</sup> TTS in 1.5 T



# **Collection efficiency in 1.5 T**



The collection efficiency decreases by ~5% in 1.5 T. The RMS is consistent with the statistical error of each measurement.

# Lifetime

## Life test setup

25

 $\Delta QE(Q)$ 


Load the output charge of the MCP-PMTs by the LED. – The output charge is measured by a CAMAC ADC. Monitor the hit rate ( $\infty$ QE) by the laser single photons.




# QE vs. output charge

QE declines as a quadratic function of total output charge:

 $\frac{\text{QE}(Q)}{\text{QE}(0)} = 1 - \left(\frac{0.447}{\text{Life}}Q\right)^2$  Life: Output charge at QE(Q)/QE(0) = 0.8








### **Prospects of performance improvement**

#### Uniform gain?

28

- Production process of MCPs may affect the uniformity.
- Better timing resolution?
  - Shorter distance from the photocathode to MCP1 would suppress the TDC tail due to the electron recoil on the MCP1 surface.





Photocathode

MCP2

- Lifetime improvement for the coated MCP-PMT.
  - R&D is ongoing for future production.

Improvement of the efficiency would be most significant in terms of the TOP PID performance.



We succeeded in mass production of
 512+α MCP-PMTs for the TOP counter.
 Achieved good quality as
 virtuous cycle of [production ->
 performance check -> production] created.

**Higher Quality** 

Systematic measurement of a large number of
MCP-PMTs for the first time:
→ Understand the performance of the MCP-PMT well.
→ Improve the performance for future production.

<u>600</u>

Performance check