

# Bead-Pull RF Measurement System


Jackline C. Koech

University of Massachusetts

Supervisor: David Peterson

Antiproton Source, Accelerator  
Division

08/09/2011



# Outline

- Bead-Pull RF Measurement System
- Motivation
- Project Setup
- LabVIEW Program
- Test Measurements


# Bead-Pull RF Measurement System

- Bead (dielectric/metallic) being pulled through a cavity
- Take Resonant frequency measurements
- Goal: Know the electric field distribution inside the cavity.
- Slater Perturbation Theory


$$\frac{\omega^2 - \omega_0^2}{\omega_0^2} = k \int_{\Delta r} \frac{\mu H^2 - \epsilon E^2}{2U} dv$$

# Motivation

- Most accelerators use RF cavities
- After fabrication, use Bead-Pull System to ensure that the electric field distribution is as designed



# Bead-Pull RF Measurement System



- Cavity
- Bead, Thread & Pulley System
- Motor
- Network Analyzer
- Computer (control software)





# Equipment

## Step Motor



- HT23-601 by Applied Motions
- 20,000 steps/rotation
- Approx. 255 steps moves the bead by 1mm


## Network Analyzer



- Agilent 8720ES
- Takes the Radio Frequency measurements

# Objectives


- LabVIEW program that controls both the Step Motor and Network Analyzer
- Check the accuracy of the measurements from the Network Analyzer and the motion of the Step Motor






- National Instruments' Software
- Graphical Programming Environment

**A LabVIEW program is called a virtual instrument (VI)**



Front Panel



Block Diagram

# Approach

## Network Analyzer VI

Control how measurements are taken

## Step Motor VI

Control the motion of the bead



## Overall Bead Pull System VI

Operates both the Step Motor VI and Network Analyzer VI as required



# BEAD-PULL RF MEASUREMENT SYSTEM

## STEP MOTOR PARAMETERS

Velocity: 10 Acceleration Rate: 25 Deceleration Rate: 25

## INPUT SPECIFICATIONS

### ANALYZER MODEL SELECTION

AGILENT 8720ES

### OUTPUT SELECTION

RESONANT FREQ MEASUREMENTS

### File Datapath

Z:\Jackline Koech\LabVIEW\HP4396\TEST.csv

Steps/mm: 254.71 Step size(mm): 0.1 Start distance: 0 Stop distance(mm): 200

### Comments

Agilent 8720ES, resonant freq measurements

NWA VISA resource in: COM4

NWA VISA resource out: COM4

Motor VISA resource: COM1

Motor Visa resource out: COM1


error in

error out

status: code: 0

status: code: 1073676294

## RESONANCE FREQUENCY MEASUREMENTS



## OUTPUT PANEL

STEPS/stepsize Distance moved Total steps moved

25 59.2 14800

Running



### Data

```
0,-4.884609E+01, 4.538155E+04, 1.487030000E+09,
1,-4.881508E+01, 4.582165E-24, 1.486970000E+09,
2,-4.878693E+01, 1.161769E+07, 1.487060000E+09,
3,-4.883217E+01, 6.386886E+18, 1.487040000E+09,
4,-4.883993E+01, 8.535699E-33, 1.487080000E+09,
5,-4.876726E+01, 6.006204E-19, 1.487050000E+09,
6,-4.883616E+01, 4.127246E-08, 1.486990000E+09,
7,-4.884594E+01, 3.630524E+05, 1.487040000E+09,
8,-4.886055E+01, 7.872094E-14, 1.486990000E+09,
9,-4.882817E+01, 8.535700E-33, 1.487130000E+09,
10,-4.882404E+01, 1.056668E-05, 1.487100000E+09,
11,-4.881540E+01, 3.579975E-26, 1.487030000E+09,
```

# Measuring the Resonant Frequency




# Pill Box Test Cavity Measurements



The bead-pull system moves the bead accurately as required

The system acquires the data correctly from the Network Analyzer



Simulated electric field

# Conclusion

We have successfully developed a software that controls the Bead-Pull System hardware

- Moves the bead accurately
- Takes Measurements accurately



# Acknowledgements

- Fermilab SIST committee
- Supervisor: David Peterson
- Ding Sun and His group
- Mentors: Cosmore Sylvester and Mayling Wong
- Dr James Davenport.
- Antiproton Source Department

# Questions??

