
ARAPUCA

Ana Machado



Work principle



The simplest geometry is a **box** with highly reflective internal surfaces with an open side.

The open side hosts the **dichroic filter** it has the property of being **highly transparent** for wavelengths **below a cutoff** and **highly reflective above it.**


The filter is deposited with **TWO SHIFTERS** – **one on**

each side


The shifter on the **external side**, S1, converts LAr scintillation light to a wavelength L1, with **L1 < cutoff**


The shifter on the **internal side**, S2, converts S1 shifted photons to a wavelength L2, with L2 > cutoff

The internal surface of the ARAPUCA is observed by one or more SiPM





Dichroic Filter





TPB and pTP emission spectrum






ARAPUCA

- Large area photon collector
- Modular and thin device easily adaptable to different mechanical constraints
- Can detect light on both sides
- Can be detect VUV or/and visible light
- Efficiency can be tuned by varying the photo cathode coverage. Currently at few percent level (3-4%) with a limited coverage with active SIPM
- Other applications

LARGE AREA PHOTON COLLECTOR

Gain factor between 7 and 10



In units of 8" PMT: 1 *X-ARAPUCA bar* is equivalent to 2.5 8" PMTs -> 25 equivalent PMTs/APA

