

How standards are *accelerating* the adoption of machine source radiation sterilization

Emily Craven, Director Global Sterility Assurance
Boston Scientific
Medical Device Sterilization Workshop
Sept 24, 2021

Identifying the need

- For the medical device industry to continue to grow, we need to provide more options for radiation sterilization
- In order to take advantage of machine source radiation, guidance can help organizations understand fundamentals of:
 - Designing products for radiation
 - Transferring products from one radiation source to another
- The collaboration between industry and regulatory and between standards organizations is helping to fill this gap

Relevant Standards Organizations


- AAMI ST – Sterilization Standards Committee
 - WG2 Radiation Sterilization
 - WG8 Microbiological Methods
 - WG15 Assurance of Sterility
 - WG96 Compatibility of Materials
- ASTM E61 Radiation Processing
 - E61.01 Dosimetry
 - E61.02 Dosimetry Systems
 - E61.03 Dosimetry Applications
 - E61.04 Specialty Applications
- ISO TC198 Sterilization of health care products



How standards help

- Recent work in standards is helping to lay the groundwork for easing the transfer between sterilization modalities
 - Ready for publication: **TIR100, *End-to-end microbiological quality and sterility assurance***
 - Ballot passed, preparing final draft for review: **TIR104, *Guidance on transferring health care products between radiation sterilization sources***
 - Work started on revision to **ISO 11137-1, *Requirements for development, validation and routine control of a sterilization process for medical devices***
 - New revision initiated: **TIR17, *Compatibility of materials subject to sterilization***
 - New ASTM standards on modality specific OQ

How standards for radiation sterilization work



TIR100 – End-to-end

- What is meant by end-to-end?



- Decisions made during all stages of product design and lifecycle can have implications for sterilization options

How does TIR100 help?

- TIR100 is not a radiation specific document but provides guidance on decision making at each of the stages that can affect the ability of a product to be radiation sterilized.



- Can I choose radiation compatible materials?

SOURCE



- Where is my sterilization capacity?

DELIVER



- How is my product quality maintained?

PLAN



- What will my sources of bioburden be?

MAKE



- What controls are in my work environment?



- Are my requirements being met?

Design for Sterilization - Examples

Materials selection and intended function

- Alternate materials, e.g. PCTFE or PVF vs PTFE or Radiation stabilized PP vs PP
- Leachables/extractables post-sterilization (plasticizers, fillers, additives, antioxidants) affecting biocompatibility
- Potential for beneficial sterilization induced changes (i.e. annealing, reducing solvents, curing hydrophilic coatings)
- Residuals as a function of mode of patient contact and patient population
- Does the product need to be sterile? (accessories, cables)

Design and Manufacturing

- Tight interferences between mated surfaces (e.g. stoppers, metal-metal, etc.)
- Areas of high density
- Bioburden controls requirements for manufacturing environment and incoming materials

Testing protocols

- Extent of verification testing allows re-sterilization and/or response to deviations in process

Packaging Considerations

- Gas permeability may impact shelf life for combination products that have oxygen sensitivity – added packaging step vs qualifying radiation process
- Orientation of product within packaging, repeatability and rigidity for radiation processes
- Amount of material going through sterilizer and impact on process efficiency and fugitive emissions
- Size of product box relative to sterilizer

Sterilization Site

- Capacity availability and back up
- Turn around time and inventory considerations
- Carbon footprint and/or extra packaging associated with transportation
- Cold chain or special environmental requirements

Regulatory efficiencies

- Time to market – novel vs established, predicate products etc.
- Product families – leveraging existing product testing data

TIR17 – Compatibility of materials

- TIR17 is a resource that can be used in new product design when selecting materials for a health care product
- Valuable update provided in 2017, but as materials evolve and new sterilization methods develop, new information should be incorporated
- Work on the next revision is underway
 - Better guidance on ranking systems for material compatibility
 - Modality selection process guidance
 - Guidance on appropriate challenge conditions and testing
 - Guidance on evaluating changing modalities?

TIR104 – Transfer between radiation sources

- TIR104 was initiated to clarify (and correct) information in ISO 11137-1 on transfer between radiation sources
- Evaluation of process capability when making the decision to transfer
 - Can the irradiator physically deliver the dose required?
 - Will dose specifications need to be updated (and what is the likelihood of success?)
- Guidance on transfer of both:
 - Sterilization dose
 - Maximum acceptable dose

Transfer fundamentals

Minimum Dose is Dose

- When product is dry, sterilization and verification doses do NOT need to be re-established
- When product is wet or can support microbial growth, a dose audit can be used to verify that sterilization dose is appropriate
 - Also applies to changing locations with same modality

Maximum Dose is not always Dose


- Dose rate and temperature may have an impact on maximum dose suitability
- In general transfers from low dose rate to high will not require retesting unless a new maximum dose needs to be established to meet DUR requirements
- Activation assessment required for e-beam $>10\text{MeV}$ or x-ray $>5\text{MeV}$

TIR104 – End-to-end considerations

- How easy it is to transfer dose depends on how much information you already know about your product
 - Sterilization dose: How did you determine? What bioburden controls do you have in place to reduce minimum dose?
 - Maximum acceptable dose: Do you know your REAL maximum dose or did you just challenge the process you had available?
- Does your product rely on the sterilization process for a functionality enhancement vs failure?
 - Crosslinking that makes product or packaging stronger?
 - Heating which anneals or cures

Transfer – How design elements help

- Packaging: Does your packaging design allow for a presentation to both electrons and photons?
 - Try not to have high density areas overlap
 - Try to keep nuclear thickness as uniform as possible
 - Consider both single unit and shipper configurations
 - Will your package size work at multiple irradiators (i.e. pallet vs tote vs conveyor)
- Is your packaging and product together designed to allow for heat transfer out of the package?
- Can your product be flipped for a vertical electron beam process?



Other work in standards

ISO 11137-1 revision


- Align with ASTM standards and ISO 11137-3 and -4
- Align with ISO 11137-2 and 13004
- Provide more specific allowances and guidance for parametric release
- Look at changing the threshold for activation assessment
- Align with TIR104

ASTM new documents

- Provide better technology specific OQ guidance
 - Gamma completed
 - Electron beam in process
 - X-ray next!
- Update and author standards as appropriate to be referenced in revised ISO 11137-1

Other publications

- Peer reviewed publications can lay the groundwork for future guidance and standards
- Two recent groups of AAMI BI&T supplements organized by the Kilmer Collaboration team, article are referenced in TIR104:
 - Industrial Sterilization, Process Optimization and Modality Changes: www.aami.org/sterilization-supplement-2020
 - Industrial Sterilization, Challenging the Status Quo, Driving for Continuous Improvement: <https://www.aami.org/news-resources/publications/bi-t/is-supplement/is-supplement-2021>
- Articles published in Radiation Physics & Chemistry on specific topics also used to support TIR17, ISO11137-1 including outputs from Team Nablo and IMRP



In summary

- New guidance will help **accelerate** the adoption of machine source radiation sterilization technology
- Consideration of **end-to-end** product lifecycle will make modality transitions easier
- **Coordination** between different standards organizations is key
- Standards work is a great example of **collaboration** to meet the needs of the health care product sterilization community
 - There are opportunities to **get involved**, please take advantage!