

Upgrades and Partnerships Enhancing High-Dose Radiation Calibration at NIST

Today's topics

Calibration
services

Lab renovation
and move

Acquisition of
new irradiator

Consensus value
to scale e-beam


ERDS
collaboration

NOAC
collaboration


NIST High Dose calibration services



- Calibration Irradiations of Customer Supplied Dosimeters with Co-60 Gamma Rays (4901XC)



- Dose Interpretation of Customer-Irradiated NIST Transfer Dosimeters (4902XC and 4903XC)



Dosimetry for High-Dose Applications

High-Energy Electron Beams

Photon Beams

25 Results



Sequence Name Newest Price



Irradiation of a customer supplied dosimeter with Co-60 gamma-rays

SKU: 49010C

Unit of Measure: Each

Price: \$2,328.00

Product Search

Browse NIST products by name.

Search Products



My Cart

You have no items in your shopping cart.

Irradiation of a customer supplied dosimeter with Co-60 gamma-rays



Irradiation of a customer supplied dosimeter with Co-60 gamma-rays

SKU: 49010C

Availability:

NIST has temporarily suspended this service. To inquire about placing an order, or for more information, please contact the primary technical contact listed below.

Primary NIST Technical Contact:

Name: Ileana Pazos

Phone: (301) 975-4121

Email: [Email NIST Technical Contact](#)

Secondary NIST Technical Contact:

Name: Lonnie Cumberland

Phone: (301) 975-6869

Email: [Email NIST Technical Contact](#)

Price: \$2,328.00


On March 27, 2019, NIST introduced a new online shopping experience to order Calibration services. Click here to access the new [NIST storefront](#)

Renovated EPR facility




- 4 EPR spectrometers and support instruments moved to renovated space.
- Services are suspended during the move, setup and testing.
- Revalidation of the service is in progress.

EMX enhancements




÷




÷



EPR calibration curves



Automatic pellet exchange system



Virtual summer student project: vacuum pen, pellet tray, pellets, and a robot arm were sent to the student's home.

Next steps: setup near EMX and integrate in auto-collection.

Dose Rate Effect

Limiting upper dose for alanine dosimeters

Rate effect is estimated to be

- zero at dose rates above 2 Gy/s
- significant at some value below 2 Gy/s, and
- clearly measurable at 1 Gy/s

Rate effect depends on absorbed dose, it

- is not measurable at 1 kGy or less
- becomes significant above 5 kGy
- reaches a maximum effect at about 30 kGy

The effect is intrinsic to alanine and is not dependent on the chemical form or manufacturing formulation of the alanine dosimeter.

The study postulated that the production of one (or more) of the radiation-induced alanine radicals is dependent on the dose rate.

[Desrosiers et al., 2008](#) and [Desrosiers and Puhl, 2009](#)

<https://www.nist.gov/programs-projects/basic-metrology-dosedose-rate-effects-alanine-dosimetry>






Fig. 1. A quality-control check-standard plot for the high-dose transfer dosimetry service. The relative difference of the computed absorbed dose for the check-standard measurement from the absorbed dose delivered is plotted versus the measurement date. Three doses were measured on each date: 1 kGy, 10 kGy, and 50 kGy (the symbols for each are defined in the legend). Dosimeters that measured $\pm 1\%$ from the target dose are deemed acceptable.




Acquisition of new irradiator



The Hopewell Designs, Inc. Model GR420 High Dose Rate Self-Contained Research Irradiator is a modern replacement for the legacy systems which are no longer supported.

The GR420 combines the well-known and understood capabilities and features of legacy systems with numerous state-of-the-art design upgrades including larger chamber volume, a rotating chamber, higher precision timing/process control and greatly improved external (operator) dose rates.

- NIST high-dose calibrations and irradiation services are performed in GammaCell Co-60 irradiators (designated **GC**).
- Metrology is based on alanine-EPR dosimetry, traceable to the primary standard water calorimeter under the Co-60 vertical beam (**VB**).




Comparisons



CCRI-S3

The Consultative Committee for Ionizing Radiation (CCRI) third supplementary comparison (S3) of standards for absorbed dose to water in Co-60 gamma radiation was completed

7 laboratories participated



New NIST irradiator

Dose rate determined by comparison to NIST irradiators (< 5 kGy)
International comparison for higher dose will be completed

validate GR420



ASTM Proficiency Test

Proficiency Test Plan for the Characterization and Evaluation of High-Dose Dosimetry Techniques in Radiation Processing

industry comparison

Revalidation plan



Phase I: March –TBD

Suspend HD T&G
(all doses)

**Move EPR lab to
renovated space;
Set up
spectrometers**

Update Hazard
Reviews/SOPs

Re-establish
calibration curves
< 5 kGy
Using legacy
irradiators

Revalidate
restart T&G < 5 kGy

Phase II: TBD

**Irradiator lab # 2
renovation complete
and receive GR420**

Update Hazard
Reviews/SOPs

Initial safety and
uniformity testing

Establish dose rate &
traceability. Determine
geometry &
temperature
correction factors;
Comparisons

Validate and start
T&G < 100 kGy

Phase III: TBD

**Irradiator lab # 1
renovation begins**

Partial operation

✓ T full range

✓ G > 1 kGy *

✗ G < 1 kGy *


Update Hazard
Reviews/SOPs

HD - high dose
T - transfer calibration service
G - irradiation calibration service
* - lower dose range TBD

Phase IV: TBD

**Renovation of
irradiator lab # 1 is
completed**

Fully operational



EPR Lab



Irradiator Lab #1



Irradiator Lab #2

Scaling factor for high energy electron beam

Review of literature from the last 15 years on scaling from Co-60 to electron beams.

Recommendation by the lead laboratories providing alanine calibration services.

An overall conversion factor of 1.014 (standard uncertainty = 0.5 %) was obtained.

Alanine dose readings determined from a Co-60 alanine calibration should be multiplied by this factor to yield the dose delivered in an electron beam.



Radiation Physics and Chemistry


Volume 171, June 2020, 108673



Determination of a consensus scaling factor to convert a Co-60-based alanine dose reading to yield the dose delivered in a high energy electron beam

Malcolm McEwen ^a  , Arne Miller ^b, Ileana Pazos ^c, Peter Sharpe ^d

[Show more](#) 

 Add to Mendeley  Share  Cite

<https://doi.org/10.1016/j.radphyschem.2019.108673>

[Get rights and content](#)

2021 CIRMS meeting

Industrial Applications Subcommittee

(chair: Ileana Pazos, PhD, National Institute of Standards and Technology)

Current Solutions and Needs in Dosimetry for Radiation Processing Applications

Mike Pageau, GEX Corporation


Chip-Scale, Integrated Photonics for Radiation Dosimetry and Calorimetry

Nikolai Klimov, PhD, National Institute of Standards and Technology (NIST) Gaithersburg, MD

Recent Advances Towards a Deployable Emergency Response Dosimetry System (ERDS)

Jason Campbell, PhD, National Institute of Standards and Technology (NIST) Gaithersburg, MD

Emergency Radiation Dosimetry System




- NIST and Global Resonance Technologies LLC are developing a robust, sensitive and portable device to rapidly measure the amount of radiation received by exposed individuals in the aftermath of a mass-exposure nuclear event.
- Small embedded dosimeters in ID cards are a cost-effective way of distributing large numbers of dosimeters as a pre-event mitigation strategy.




- Smaller alanine dosimeters coupled with an EPR system that can detect 1 – 5 Gy could be adapted to industrial needs.
- Plastic encapsulation could be coded and tracked.
- The new dosimetry system has been demonstrated to be suitable for electron beam irradiations.

NIST on a Chip



- By measuring the temperature rise in silicon devices due to radiation absorption, NIST researchers have taken the first step toward a miniaturized version of the conventional calorimetry technique.
- The chip-based device detects how heat energy affects the properties of light passing through microscopic channels being irradiated by gamma rays or an electron beam.

<https://www.nist.gov/noac/technology/radiation/photonic-sensors>



- By leveraging nanomanufacturing and telecommunications technology, the new sensor promises highly accurate readings at smaller dimensions.
- Exploring different doping strategies that would result in accumulated defects would convert this calorimeter into a dosimeter.

STAY IN TOUCH

CONTACT US



NIST.gov



@nist