
Tolerancing Tera-BREAD Using Ray-Tracing Simulations

Matthew Malaker, FNAL for the BREAD Collaboration

FERMILAB-POSTER-21-106-STUDENT

What is BREAD?

BREAD is a dish-antenna style axion detection experiment currently in the design stage. It is a concept that allows broadband axion detection across a wide range of frequencies [1]. Axions interact with magnetic fields in the presence of conductors, emitting photons normal to the event surface in the process. BREAD is designed to fit into a solenoid magnet and focus the axion light to a focal spot. We currently plan on using an SNSPD, which has a dark count rate of around 1/day.



Concept art of BREAD. The instrument consists of a revolved parabolic reflector in a cylindrical container. BREAD is a dish antenna wrapped into a cylinder. We designed BREAD this way in order to use readily available solenoid magnets. We plan on a 20cm radius reflector. Fig. Courtesy Stefan Knirck (FNAL)





Image of a superconducting nanowire single photon detector (SNSPD), which are typically 400μm*400 μm [2].


[Hochberg et al., PRL 123 (2019)]

This manuscript has been authored by Fermi Research Alliance, LLC under Contract No. DE-AC02-07CH11359 with the U.S. Department of Energy, Office of Science, Office of High Energy Physics

they do so.


Studies of Surface Ripples Using FRED





(Left) BREAD reflector with exaggerated ripple in Z direction. (Right) Reflector Z Ripple focal spot. As perturbations increase, these spots increase in size, but their features scale uniformly. Irradiance is multiplied by the radius to extract more information, such as divergence at the origin.

Sinusoidal perturbations (ripples) cause a broadening of the focal spot. The width of the focal spot is determined by the product of the ripple angular frequency and amplitude. To first order, degeneracies in this parameter space yield the same spot size. This means that these two surface quality measures are coupled in terms of tolerances and must be treated accordingly.



Each curve is a set of surface quality parameters that give a spot size of 0.5mm in radius. Each curve is a different ripple, either in a different direction or on a different surface (e.g emitter or reflector)

We found that BREAD is significantly more sensitive to angular perturbations than positional ones. FRED defines surface positions and reflection angles separately, which enabled us to test the outcome of removing the position component of a perturbation. As such, the tolerances of the emitter spatial dimensions are much looser than any ripple. Surface ripples can cause significantly larger deflection angles than other perturbations and must be more tightly constrained as a result.

| Perturbation | Formula | Example Parameters |
|--------------------|--------------------------------|--|
| | A | w _ A 363mm |
| Z-ripple (Emitter) | $w = 9079 * \frac{1}{\Lambda}$ | $\overline{0.5mm} - \overline{20\mu m}^* \overline{\Lambda}$ |
| Z-ripple | W - 7690 r | w = A * 307 mm |
| (Reflector) | $w = 7680 * \frac{1}{\Lambda}$ | $0.5mm - 20\mu m$ Λ |
| Φ-Ripple | w - 220000 # | w = A 9.2 Deg |
| (Emitter) | $w = 230000 * - \varphi$ | $0.5mm - 20\mu m^{*} \varphi$ |
| Φ-Ripple | $w = 190000 * \frac{A}{-}$ | w A _ 7.4 Deg |
| (Reflector) | $\psi = 190000 * \varphi$ | $0.5mm - 20\mu m^{\uparrow} \varphi$ |
| Uniform Radial | 1 00 . [47] | $\frac{w}{} = \frac{R}{}$ |
| Deviation | $w = 1.98 * \Delta R $ | 0.5mm 	 0.25mm |
| Conical Radial | $w = 3.39 * \Delta R $ | <u>w</u> R |
| Deviation (Top) | | 0.5mm = 0.15mm |
| Conical Radial | | w R |
| Deviation | $w = 1.39 * \Delta R $ | $\frac{1}{0.5mm} = \frac{1}{0.36mm}$ |
| (Bottom) | | |

Conclusions

In conclusion, we determined the tolerances to which the Tera-BREAD instrument must be fabricated. This includes the basic dimensions of each piece as well as constraining macroscopic ripples on each surface. Such ripples take the form of a sinusoidal perturbation, which showed that this design is more sensitive to deflection angles than to those in spatial coordinates.

References

- [1] Horns et al, JCAP, Issue 04, article 016. arXiv:1212.2970
- [2] Hochberg et al., PRL 123 (2019)
- [3] FRED: https://photonengr.com/fred-software/





In order to deduce mechanical

tolerances, we used FRED, a

Mechanical imperfections will

defocus the beam, which must

be kept in check because of the

small size of SNSPDs, pictured

on the left. This feature is the

reason for this study. Our goal

mechanisms which broaden the

spot size and in which ways

was to determine the

raytracing software package [3].