Cold Neutron Moderator Testing at LENS

D. V. Baxter, W. M. Snow, C-Y Liu Indiana University/CEEM










CEEM: Areal View



CEEM: Cutout View/Facilities




Proton Radiation Therapy (>1000 patients so far)

Low Energy Neutron Source (LENS)

- Based on **proton linac** and low-energy (p,nx) reactions (E_p <13MeV) in Be. 2 target stations:
- Unmoderated target for MeV neutron applications in science and engineering (Single Event Effects in electronics, defense, radiography).
- Flexible **target-moderator-reflector** system for moderator development (solid CH₄ at 4K)
- Variable pulse width (from $<5 \mu s$ to more than 1.0 ms).
- In "long-pulse" mode, LENS has time-averaged cold neutron intensity suitable for **neutron scattering** (SANS, neutron spin echo).

LENS Missions (our Venn Diagram)




LENS Accelerator

~20 mA peak proton current, Variable frequency/pulse width


13MeV proton linear accelerator (RFQ/DTLs)

High-power Klystrons in operation

"Neutron physics on a human scale" (F. Mezei)




LENS Facility Layout




Designed/built/characterized by graduate students Local user program in operation

Target Moderator Reflector (TMR)




Target/Moderator Test Assembly




PT-410 pulse tube cryorefrigerator

Top view of target/moderator



"extra" space for flexibility

4K LENS Cold Neutron Moderator




Research on Prototype Moderators






View of the reflector (inside a lead cask to shield gammas) and the cavity available for test moderators. The proton beam enters from bottom of left-hand image.

Fit to the Spectrum 13 MeV, CH4, 6K



Phase II Structure of Solid Methane





CH4 is the brightest-known cold neutron moderator (not useable at high power sources due to rad damage)

Below 20 K 1/4 of the CH4 are free rotors 3/4 are hindered rotors

Above 20K All sites are free rotors

S(q,omega) in CH4 Phase I and II (Y. Shin)







Phase II


Phase I

PE/Si vaned moderator test




Test idea (supported by MCNP) of improved extraction of cold neutrons from deep within the moderator volume (Baxter et al, Physics Procedia [Elsevir, 2012])







Diamond nanoparticles as reflector?

700







R Cubitt et al. NIMA 622, 182-85 (2010)

Can high-albedo materials such as diamond nanoparticles be used as radiation-hard reflectors near the moderator to improve cold/VCN brightness?


Magnetic Shielding for nnbar

W. M. Snow Indiana University/CEEM





N-Nbar search at ILL (Heidelberg-ILL-Padova-Pavia)



Baldo-Ceolin M. et al., Z. Phys. C63,409 (1994).

Quasifree Condition: B Shielding





Fig. 10. The transverse field compensation system. Loops 1 and 2 are under 49 A current and compensate the horizontal field component; loops 3 and 4 are under 120 A current and compensate the vertical field component.

Magnetic Shielding Measurements for neutron spin rotation measurements at NIST





B~1 nT are routinely obtained in relatively small (few m³) volumes with multiple shields

Large size of any n-nbar apparatus poses a challenge.

Magnetic Shielding for NEDM Experiment At Petersburg Nuclear Physics Institute



U. Schmidt et al. / Propagation of a neutron beam

Quasifree Condition: Neutrons as the magnetometer

Performed by polarized neutrons and use of neutron spin echo spectroscopy (U. Schmidt et al, NIM A320, 569 (1992).

May need polarizers and analyzers with larger phase space acceptance to polarize and analyze beam






Fig. 3. Polarized neutron analyser equipment at the exit of the zero-field region, at 70 m distance to the polarizers of fig. 2.





Fig. 4. Neutron-spin echo curve produced with the apparatus of figs. 2 and 3. The zero-offset of this curve measures the residual magnetic field along the magnetically shielded 70 m long neutron beam line to $\langle B_z \rangle = 4$ nT, via a neutron spin precession angle of $\phi = 2^\circ$. The solid line is a fit to the theoretical lineshape.

POLARIZED ³He for Neutron Polarimetry

Use neutrons as magnetometers. Polarize/analyze neutron beam using 3He



Large neutron phase space acceptance Polarizer/analyzer pair can measure B using neutron spin rotation


Required R&D

- Identification of scheme for suppression of external magnetic fields/field gradients in the presence of the μ-metal;
- Identification of an assembly, annealing, and demagnetization strategy;
- Measurement of the ambient magnetic field at the proposed location;
- Preliminary mechanical design of the shield;
- Analysis of the loss of sensitivity of the oscillation measurement to residual magnetic fields;
- Investigation of materials compatibility for the vacuum chamber.
- Verify maintenance of the field condition


Supermirrors (multilayer mirrors)

H. Shimizu

Experimental Setup of ILL nn Experiment



Neutron Optics for nn Oscillation



Neutron Supermirrors with m=3.5 are "Items of Commerce"




Focusing Guide



Guide Array in Reactor




Guide Installation





Guide Manufacturers Catalog

Photos Swiss Neutronics, LANL, ILL










critical angle of total reflection


φ_c(Ni)/λ=1.7 mrad/Å v_⊥(Ni)=7 m/s

Multilayer Mirror (Monochromator)




φ_c(Ni)/λ=1.7 mrad/Å v⊥(Ni)=7 m/s

Multilayer Mirror (Monochromator)



φ_c(Ni)/λ=1.7 mrad/Å v⊥(Ni)=7 m/s


Supermirror



 $m = \phi_c/\phi_c(Ni) = v_c(Ni)/v_c$

φ_c(Ni)/λ=1.7 mrad/Å v⊥(Ni)=7 m/s

Supermirror




 $m = \phi_c/\phi_c(Ni) = v_c(Ni)/v_c$

φ_c(Ni)/λ=1.7 mrad/Å v_⊥(Ni)=7 m/s


First m=7 Supermirror

m=5 (4766 bilayers) + m=7 (7131 bilayers)




m=4-7 Supermirrors

Supermirror: commercially available up to m=7 (v_⊥=50m/s)




Option B 10-100m Vertical Flight Path


n-nbar experiment using very cold neutrons



n-nbar experiment using very cold neutrons



Folding acceptance and spectrum: illustration



(Some) Issues

- Current neutron guide coatings are optimized for many reflections: High reflectivity more important than high critical angle.
- Current guide construction methods are not easily "scalable"
- Current world wide capacity is ~1000 m²/year
- No commercial source in the US.

•...

Summary

Multilayer mirrors can enhance the figure-of-merit for n-nbar experiments.

Multilayer fabrication technology was remarkably improved in the past decade.

supermirrors m≤7 exist

Other (coupled) issues:

Focusing geometry (~elliptical)
Spectrum