# Advantages of Operating your TPC at the Drift Velocity Maximum and the Properties of P2


Philip Hamacher-Baumann

ND-GAr: HPgTPC+ECAL Weekly Meeting

2022



#### The DUNE Near Detector



#### ND-LAr

- Liquid Argon TPC
- Not magnetized

#### ND-GAr

- High pressure TPC
- Baseline mixture P10 (Ar + 10% CH<sub>4</sub>)
- Magnetized + contains an ECAL
- Spectrometer for particles leaving ND-LAr
- Additional v-Ar measurements
- Very low tracking threshold


#### **SAND**

- Permanently on-axis beam monitor
- Magnetized






## **Properties of P10 under ND-GAr Conditions**



Maximum in  $v_d$  at 130 V/cm/bar = 1300 V/cm at 10 bar



Diffusion very low across all fields





#### **TPC Gas in DUNE's ND-GAr Detector**



Maximum in  $v_d$  for most stable for operation, at 2.5 m drift:

P10: 130 V/cm/bar -> 325 kV on cathode



P2: 50 V/cm/bar -> 125 kV on cathode



And twice that for single anode 5m drift.

Drift distance 5m or 2.5m





#### **TPC Gas in DUNE's ND-GAr Detector**



Maximum in  $v_d$  for most stable for operation, at 2.5 m drift:

P10: 130 V/cm/bar -> 325 kV on cathode



P2: 50 V/cm/bar -> 125 kV on cathode



And twice that for single anode 5m drift.

Maybe we don't need to operate at the drift velocity maximum then?

Drift distance 5m or 2.5m





## Advantages of Operating at the Drift Velocity Maximum





## **Short Reminder: Density Scaling in Gaseous Detectors**

## Most experienced at atmospheric conditions

Swarm parameters depend on T and p (+ also E)

- averages of microscopic e<sup>-</sup> atom/molecule collisions
- change in mean free path affects all parameters

All swarm parameters depend on the gas density.

magnitude	scaling $(n = N/N_0)$
electron, ion drift velocity $v_d$	$v_d(E/n)$
electron, ion diffusion coefficients $D_{L,T}^*$	$\frac{1}{\sqrt{n}}D_{L,T}^*(E/n)$
attachment coefficient $\eta$	$n \cdot \eta(E/n)^{*a}$
Light transparency $\mathcal{T}$	$\exp\left(-n\Pi_a L^*\right)$
scintillation probability $P_{scin}$	$\frac{1}{1+n auk}$
particle range $R$	R/n
Fano factor $F_e$ , $W_I$ , $W_{ex}$	$\sim { m constant}$
charge multiplication coefficient $\alpha$	$n \cdot \alpha(E/n)^{*b}$
secondary scintillation coefficient $Y$	$n \cdot Y(E/n)^{*b}$

arxiv:1710.01018[3]







#### Pressure control:

No problem – TPC in pressure vessel



High flow, high precision dome regulators



Various spring-loaded regulators





#### Pressure control:


No problem – TPC in pressure vessel

#### Temperature control:

Electronics / ECAL *inside* of vessel Could mean large temperature gradients

#### Q:

Without any temperature control, what is the spread of  $v_d$  at a constant E/p? Assume T from 0°C to 60°C







#### Pressure control:

No problem – TPC in pressure vessel

#### Temperature control:


Electronics / ECAL *inside* of vessel Could mean large temperature gradients

#### Q:

Without any temperature control, what is the spread of v<sub>d</sub> at a constant E/p? Assume T from 0°C to 60°C

#### A:

Basically negligible at v<sub>d</sub><sup>max</sup>! Ignoring T elsewhere can add *O*(10cm) error!



Height of v<sub>d</sub> curve at constant E/p is resulting temperature error.





#### Pressure control:

No problem – TPC in pressure vessel

#### Temperature control:


Electronics / ECAL *inside* of vessel Could mean large temperature gradients

#### Q:

Without any temperature control, what is the spread of  $v_d$  at a constant E/p? Assume T from 0°C to 60°C

#### A:

Basically negligible at v<sub>d</sub><sup>max</sup>! Ignoring T elsewhere can add *O*(10cm) error!



Height of v<sub>d</sub> curve at constant E/p is resulting temperature error.





## Perks of Operation at Drift Velocity Maximum: Impact of Outgassing

1.5 bar P10


Estimating effect of gas purity in HPGMC:

- Completely shut off gas
- Repeat measurements
- Watch for deterioration

Outgassing accumulates quicker at low pressures

- At 1.5bar, the change in  $v_d$  is:
  - -20% at lowest field
  - -2% at v<sub>d</sub>max
  - +2% for highest field

Sum of e-gas xsec *lowest* around v<sub>d</sub><sup>max</sup>, yet influence of impurities is **not higher**, but lower through this!







## Perks of Operation at Drift Velocity Maximum: Impact of Outgassing

1.5 bar P10


Estimating effect of gas purity in HPGMC:

- Completely shut off gas
- Repeat measurements
- Watch for deterioration

Outgassing accumulates quicker at low pressures

- At 1.5bar, the change in  $v_d$  is:
  - -20% at lowest field
  - -2% at v<sub>d</sub>max
  - +2% for highest field

So how can we get ND-GAr to operate at a v<sub>d</sub> maximum?







# Ar:CH<sub>4</sub> Mixtures with low CH<sub>4</sub> Fraction





## Properties of Argon: Methane Mixtures at E<sub>drift</sub>(v<sub>d</sub><sup>max</sup>): Drift Velocity

#### Which value do we want?

Common shape of Ar:CH<sub>4</sub> mixtures

- CH<sub>4</sub> content scales v<sub>d</sub> and shifts E/p
- High  $CH_4$  = high  $v_d^{max}$
- And vice-versa


Drift velocity determines "pile-up" between spills

- Beam arrives at ~1Hz
- Spills never overlap: 1<sup>μm</sup>/<sub>ns</sub> = 1000<sup>m</sup>/<sub>s</sub>

Each spill has a substructure

- 1µs spacing (my guess?)
- Separatable with  $v_d = 2500(5000)^{\mu m}/n_s = unreachable$

No real  $v_d$  criteria derivable from beam, can go solely by cathode voltage argument (2%  $CH_4$ ).







## Properties of Argon: Methane Mixtures at $E_{drift}(v_d^{max})$

Like E/p scaling, magnetic field also scales B~B/n


- 10 bar and 0.5 T only reduced d₁ by ~1%
- B field still mandatory for sign selection

Longitudinal diffusion below ALICE reference

Transversal diffusion is larger:

- Point charge spread after 250cm (500cm):
   5.5 mm (7.8mm)
- ½ to 1 pad size


Diffusion heavily suppressed by pressure scaling - even for P2 and low fields.








## **Properties of P2 under ND-GAr Conditions**






Diffusion very low across all fields





## Gas Gain on Wire(s)

- Townsend multiplication happens a few wire-radii above wire
- ALICE Anode Wires have 20µm diameter
- MWPC and single wire look the same 100µm scale
  - Use a Single Wire Proportional Counter as proxy



magnitude	scaling $(n = N/N_0)$		
electron, ion drift velocity $v_d$	$v_d(E/n)$		
electron, ion diffusion coefficients $D_{L,T}^*$	$\frac{1}{\sqrt{n}}D_{L,T}^*(E/n)$		
attachment coefficient $\eta$	$n \cdot \eta(E/n)^{*a}$		
charge multiplication coefficient $\alpha$	$n \cdot \alpha(E/n)^{*b}$		
secondary scintillation coefficient $Y$	$n \cdot Y(E/n)^{*b}$		

arxiv:1710.01018 [3]

#### First Townsend coefficient

- Both electric field and value scaling
- 10x pressure does not mean 10x anode voltage scaling

Gas gain **not** linearly proportional to pressure.





Ye et al. "Gas amplification in high pressure proportional counters" [1]




Fig. 4. Gas gain, using the current method, as a function of applied anode voltage for argon-methane (98:2) at various pressures.




Fig. 7. Anode voltage required for a gas gain of  $10^3$  as a function of the gas pressure.

- Maximum gain appears constant up to >20 bar
- Constant gain under 10x pressure increase with ~2.5x voltage increase
- Lower-CH<sub>4</sub> mixtures require less anode voltages





P5

### **Additional Benefit**

Less methane, less inflammable!

Ar:CH<sub>4</sub> in air

- below ~ 4% CH<sub>4</sub> or
- above ~17% CH<sub>4</sub>
   will not ignite.

Fermilab has its own rules.

**Table 12.1** Limits of inflammability of some gases and vapours in air that are of interest for particle detectors, from two different sources (vol. % of total mixture)

Compound	Formula	Limits of inflammability of gas or vapour			
		[LEW 61]		[CER 96]	
		Lower	Upper	Lower	Upper
Methane	CH <sub>4</sub>	5.3	15.0	4.4	16.9
Ethane	$C_2H_6$	3.0	12.5	2.4	14.6
Propane	$C_3H_8$	2.2	9.5	1.8	10.4
Butane	$C_4H_{10}$	1.9	8.5	1.4	8.9
Isobutane	$i-C_4H_{10}$	1.8	8.4	1.55	8.4
Pentane	$C_5H_{12}$	1.4	7.8		
Isopentane	$i-C_5H_{12}$	1.4	7.6		
2,2-dimethylpropane	$C_5H_{12}$	1.4	7.5		
Hexane	$C_{6}H_{14}$	1.2	7.5		
Ethylene	$C_2H_4$	3.1	32.0		
Propylene	$C_3H_6$	2.4	10.3		
Methylalcohol	CH <sub>3</sub> OH	7.3	36.0		
Isopropylalcohol	$C_3H_7OH$	2.0	12.0		
Diethylether	$C_4H_{10}O$	1.9	48.0		

From [2]





#### **Conclusion**

Many Only benefits of operating a TPC at the drift velocity maximum!

- Temperature and pressure fluctuations almost negligible
- Diffusion low due to pressure scaling
- Resilience of swarm parameters against impurities

Reducing the CH<sub>4</sub> content in the DUNE ND-GAr gas mixture puts v<sub>d</sub><sup>max</sup> at an achievable cathode voltage

- [+] P2 has earlier onset of gas multiplication than P10
- [+] Large Ar content in both cases, dE/dx should look very similar
- [-] P2 shows no primary scintillation

#### Source(s):

- 1. Z. Ye, R. Sood, D. Sharma, R. Manchanda, and K. Fenton. Gas amplification in high pressure proportional counters. Nucl. Instrum. Meth. A, 329(1-2):140–150, May 1993. doi: 10.1016/0168-9002(93)90929-c.
- 2. W. Blum, W. Riegler, and L. Rolandi. Particle Detection with Drift Chambers. Particle Acceleration and Detection. Springer, 2nd edition, 2008. ISBN 9783540766834. doi: 10.1007/978-3-540-76684-1.
- 3. D. González-Díaz, F. Monrabal, and S. Murphy. Gaseous and dual-phase time projection chambers for imaging rare processes. Nucl. Instrum. Meth. A, 878:200–255, jan 2018. doi: 10.1016/j.nima.2017.09.024.






# Thank you!






## **Properties of 1% CF<sub>4</sub> in Ar under ND-GAr Conditions**



Maximum in  $v_d$  at 105 V/cm/bar = 1050 V/cm at 10 bar




Diffusion very low across all fields (values below 5 V/cm/bar suspicious)





## VUV Absorption Cross-sections of CH<sub>4</sub> and CF<sub>4</sub>



Absorption cross sections of methane CH<sub>4</sub> at room temperature (110-165 nm)


VUV absorption cross sections of tetrafluoromethane CF<sub>4</sub> at room temperature

From http://satellite.mpic.de/spectral\_atlas/cross\_sections/





## **Electron Scattering Cross Sections: Argon**







## **Drifting Electrons in Pure Argon**


## (Don't do it.)

Pure argon is a very slow drifter.

Electrons gain energy unhindered

Large diffusion

Excitation and ionization limit energy gain eventually.





NB: Gas multiplication in pure Ar is unstable.






## **Electron Scattering Cross Sections: Methane**





Additional rotational and vibrational degrees of freedom:


- Inelastic collisions = cooling effect
- Threshold ~ 0.2 eV








## **Electron Scattering Cross Sections: Argon + Methane**







## **Electron Scattering Cross Sections: Argon-based Mixtures**







## **Drifting Electrons in Pure CH<sub>4</sub>**


Pure methane is a very fast drifter, but only at high fields

- Saturated v<sub>d</sub> (maximum)
- Very interesting for TPC operators

Electrons cooled to thermal limit up to ~100 V/cm

Lowest diffusion

Avalanche multiplication possibly, but at very high fields.







#### **Effect of Additional Quencher**

- Electrons fall into Ramsauer minimum at higher fields
- Shift of v<sub>d</sub><sup>max</sup> to higher E
- Same behaviour could be caused by accumulating impurities with low-threshold inelastic xsec contributions
- Unclear what impurity present
- H<sub>2</sub>O and N<sub>2</sub> in sensible concentrations cannot explain shift

