

IMPLEMENTING THE CC-USB CONTROL MODULE FOR USE IN CAMAC CRATES AT THE FERMILAB TEST BEAM FACILITY

BY KAREN LIPA, SIST INTERN

AUG 8, 2012




OUTLINE

- Background info
- Specific use: cosmic ray test stand
- Goals
- My Process
- Questions



HARDWARE: CAMAC SYSTEM

Slot number: from 1 to 25



*examples on next slide



N, A, AND F

| N | A | F | Function | Data |
|--------|----|----|--|-------|
| 0 | * | 16 | Write a 16-bit marker word into the output data stream | 16 |
| 1...24 | * | * | Executes N(1..24) A(*) F(*) command on CAMAC data way | 16/24 |
| 25 | 0 | 0 | Read Firmware ID | 32 |
| 25 | 1 | 0 | Read Global Mode | 16 |
| 25 | 1 | 16 | Write Global Mode | 16 |
| 25 | 2 | 0 | Read Delays | 16 |
| 25 | 2 | 16 | Set Delays | 16 |
| 25 | 3 | 0 | Read Scaler Readout Control | 24 |
| 25 | 3 | 16 | Write Scaler Readout Control | 24 |
| 25 | 4 | 0 | Read User LED Source Selector | 32 |
| 25 | 4 | 16 | Write User LED Source Selector | 32 |
| 25 | 5 | 0 | Read User NIM Output Source Selector | 32 |
| 25 | 5 | 16 | Write User NIM Output Source Selector | 32 |
| 25 | 6 | 0 | Read Source Selector for User Devices | 32 |
| 25 | 6 | 16 | Write Source Selector for User Devices | 32 |
| 25 | 7 | 0 | Read Timing for Delay & Gate Generator A | 32 |
| 25 | 7 | 16 | Write Timing for Delay & Gate Generator A | 32 |
| 25 | 8 | 0 | Read Timing for Delay & Gate Generator B | 32 |
| 25 | 8 | 16 | Write Timing for Delay & Gate Generator B | 32 |
| 25 | 9 | 0 | Read LAM Mask | 24 |
| 25 | 9 | 16 | Write LAM Mask | 32 |
| 25 | 10 | 0 | Read CAMAC LAM (pseudo-register) | 24 |
| 25 | 11 | 0 | Read Scaler A (pseudo-register) | 32 |
| 25 | 12 | 0 | Read Scaler B (pseudo-register) | 32 |
| 25 | 13 | 0 | Read Extended Delays Register | 32 |
| 25 | 13 | 16 | Write Extended Delays Register | 32 |
| 25 | 14 | 0 | Read USB Buffering Setup Register | 32 |
| 25 | 14 | 16 | Write USB Buffering Setup Register | 32 |
| 25 | 15 | 0 | Read Broadcast Map (notepad register) | 24 |
| 26 | * | * | execute Broadcast A(*) F(*) on CAMAC dataway | 16/24 |
| 27 | ** | ** | Set Broad cast mask (3 sequential calls for 24 bit mask) | 24 |
| 28 | 8 | 29 | CAMAC Z | - |
| 28 | 9 | 29 | CAMAC C | - |
| 29 | 9 | 24 | Set CAMAC I | - |
| 29 | 9 | 26 | Clear CAMAC I | - |

- Read (F 0-7)
- Write (F 16-23)
- I: inhibit
- Z: initialize
- C: clear




HARDWARE: CC-USB

- Control module
- Transmits and receives data to/from crate modules
- Accessed by computer via USB cable
- Newest type of control module (others are obsolete)
- Not previously equipped for use at Fermilab



COSMIC RAY TEST STAND PROJECT

- Goal: to produce a reliable means of testing new detectors
- Wire chambers – x- and y-plane wires indicate location of particles
- Scintillators and PMT's omit signal when particle travels through (trigger)



HARDWARE


CAMAC crates

- Lots lying around, economical

CC-USB

- CAMAC parallel bus is now obsolete


Wire chambers,
scintillators, PMTs






HARDWARE: WIRE CHAMBERS

- x- and y-plane wires indicate locations where particles hit
- Sends a signal to card




HARDWARE: SCINTILLATORS + PMT'S

- scintillator (synthetic). A substance that glows when hit by high-energy particles or photons (www.freedictionary.com)
- Photomultiplier tubes: “multiply” result of scintillator hit (emit electrical signals)




THE SETUP

RACK Front View



RACK Top View




4 double unistruts 32" long
4 double unistruts 30" long




THE SETUP: CAMAC CRATE

- **Lecroy 3377 Time-to-Digital converter**
 - counts the time from a hit in the wire chamber until stop signal
- **Lecroy 4301: encoding and readout**
 - sends out stop signal to all the TDCs



TEST PROCEDURE

Pulse generator: produces mock STOP and hit signals




GOALS

- Software to support the use of the CC-USB control module
- Python wrapping (ease for user)
- Specifics relevant to cosmic ray project



STEPS I TOOK



PYTHON WRAPPING

All the benefits of Python

...

- Ease of use (syntax) **
- High level
- Object-oriented
- Script-based (testing)

With the power of C!

- Access to USB functions (examples)
- speed



HOW TO WRAP A C FUNCTION

STEP 1: WRAPPER FUNCTION

- Functions (PyArg_ParseTuple and Py_BuildValue) take in and return Python objects
- Within the wrapper function, call is made to C function defined elsewhere

```
static PyObject *  
spam_system(PyObject *self, PyObject *args)  
{  
    const char *command;  
    int sts;  
  
    if (!PyArg_ParseTuple(args, "s", &command))  
        return NULL;  
    sts = system(command);  
    return Py_BuildValue("i", sts);  
}
```



HOW TO WRAP A C FUNCTION

STEP 2: METHODS TABLE

- Defines the “nickname” that each function will be called by

```
static PyMethodDef SpamMethods[] = {  
    ...  
    {"system", spam_system, METH_VARARGS,  
     "Execute a shell command."},  
    ...  
    {NULL, NULL, 0, NULL}           /* Sentinel */  
};
```



HOW TO WRAP A C FUNCTION


STEP 3: INITIALIZATION FUNCTION

- Defines the name of the module to be called from python


```
PyMODINIT_FUNC  
initspam(void)  
{  
    (void) Py_InitModule("spam", SpamMethods);  
}
```




STEPS I TOOK



MY CONTRIBUTION: SOFTWARE



BIG MILESTONES



ACKNOWLEDGEMENTS

- Supervisor: Geoff Savage
- Mentors: Elliott McCrory, Jamieson Olsen
- Dianne Engram, SIST committee
- Dr. James Davenport



QUESTIONS?

