


WLS Gasses for High-Pressure Xenon Detectors



Victor M. Gehman for Dave Nygren

LIDINE2013

May 30, 2013


Fermilab. Batavia, IL

High-Pressure Gas Detectors

- Good for $0\nu\beta\beta$ and for dark matter!
- Excellent energy resolution (critical for $0\nu\beta\beta$, can also lower threshold for dark matter)
- Track imaging for high-fidelity track topology (enhanced background rejection)
- Small charge to light fluctuations (more precisely defined signal and background regions)
- Can possibly extract nuclear recoil track direction (would be a game changer for dark matter!)

Detector Schematic

Asymmetric TPC with “Separated functions”



Energy Resolution

A. Bolotnikov, B. Ramsey / Nucl. Instr. and Meth. in Phys. Res. A 396 (1997) 360–370






Fig. 5. Density dependencies of the intrinsic energy resolution (%FWHM) measured for 662 keV gamma-rays.

A Recent Result!




A Recent Result!



Tracking Too!


Real track from ^{137}Cs γ -ray – reconstructed with SiPMs



NEXT-DEMO IFIC, Valencia



What's the next step?

- Previous results were done with pure Xe (tracking plane had evaporated TPB)
- Track directionality would make a very strong case for direct detection of dark matter
- Most current experiments try for directionality by imaging the nuclear recoil track:
 - Very diffuse detectors (low target mass)
 - High energy threshold
 - Poor track image quality



A Different Approach!

- Use *columnar recombination* (CR) to extract track direction...
- Requires ionization electrons drift back through parent track:
 - Depends on angle between drift field and track direction
 - Other recombination types are independent of this angle



How to Maximize This Effect?


- Define “*Columnarity*,” $C = \frac{R}{r_0}$
- Represents the maximal difference in recombination from track angle
- In 0.05 g/cm³ xenon gas:
 - $R \equiv$ Nuclear recoil track range $\approx 2.1 \mu\text{m}$
 - $r_0 \equiv$ Onsager radius
(recombination distance) $\approx 70 \text{ nm}$ $r_0 = \frac{e^2}{\epsilon E_e}$
 - $e \equiv$ electron charge, $\epsilon \equiv$ gas dielectric constant,
 $E_e \equiv$ electron kinetic energy (usually taken as kT)
 - $C \approx 30$ in this case (would like $C > 10...$)

So What Do We Need?

- We have:
 - Short tracks (~ 70 nm)... **Don't lose electrons!**
 - Small signals... **Don't waste electrons or photons!**
- Lots of energy deposited from nuclear recoils goes into primary excitations, but...
 - excitations don't contribute to the CR signal!
 - Use the Penning Effect: convert excitons to ions with a molecular additive so that these can contribute to CR too!
- Bonus: the same molecule can cool the electrons, thus increasing the recombination probability


But Wait, There's More!

- Remember that we are detecting ionization electrons with electro-luminescence light, therefore...
- Poor photon collection efficiency means poor charge collection efficiency!!!
- We can achieve nearly 100% **coverage** if we cover the inside of the TPC with WLS plastic panels read out with PMTs (or APD's, or SiPMs, *etc.*)
- But most WLS plastic panels are not very efficient in VUV–300 nm light is pretty close to optimal though.
- Must shift 173 nm photons to 300 nm photons in the gas!




But Wait, There's More!

- Remember that we are detecting ionization electrons with electro-luminescence light, therefore...
- Poor photon collection efficiency means poor charge collection efficiency!!!
- We can achieve nearly 100% **coverage** if we cover the inside of the TPC with WLS plastic panels read out with PMTs (or APD's, or SiPMs, *etc.*)
- But most WLS plastic panels are not very efficient in VUV–300 nm light is pretty close to optimal though.
- Must shift 173 nm photons to 300 nm photons in the gas!




Two Birds With One Stone


- To extract the CR signal from a HPXe gas detector, we need two things:
 - Penning additive to convert excitations into ionizations
 - WLS that absorbs at 173 nm and fluoresces at ≈ 300 nm
- Provenance! Tri-methyl-amine (TMA) is a Penning gas known to fluoresce efficiently at 300 nm!
- also possible: Tri-*ethyl*-amine (TEA)




TMA




TEA




What Might This Look Like?



What Might This Look Like?




What Might This Look Like?




Shorter Term R&D


- The “TEA Pot”
- Measures basic response characteristics
 - Parallel-plate ionization chamber with optical sensing using 4 PMTs that look at the gap from the sides
- Will measure both light and charge as functions of density, electric field, and fraction of TMA/TEA



Shorter Term R&D



Shorter Term R&D



OSPREY: “Opportunities for Superior Performance in Rare Event Yields”

Conclusions

- This is a really unusual way to get at dark matter directionality
- Each step is quite plausible, but there are several unknowns to be addressed:
 - Penning efficiency of TMA?
 - Fluorescence efficiency of TMA in recombination?
 - Rate of ionic charge exchange?
 - Cooling rate of electrons after ionization?
- Initial simulations and R&D is underway!