

Electrospinner Upgrades for Nanofiber Production

Niko Black, SULI 2024 | Target Systems Department

FERMILAB-POSTER-24-0206-STUDENT

Supervisor: Sujit Bidhar

Background & Motivation

- Electrospinning is an inexpensive method for producing nanofibers with applications in accelerator targets, air filtration, and biomedicine
- Nanofiber diameter (varies functional properties) is determined by spinneret to substrate separation distance and polymer solution viscosity
- Current setup is incapable of adjusting spinneret displacement
- Polymer solution viscosity has not been quantified and can be applied to channel design once measured




Fig 1: Free surface electrospinner unit

Inset: Electro-spinneret channel (grooved aluminum bar)

Objectives

- Research and develop a vertical lift system for the electro-spinneret channel (positive electrode)
- Redesign the solution delivery system for enhanced reliability and precision
- Calibrate viscometer and measure the viscosities of various polymer solutions





Fig 2 (above): Scissor lift model in NX CAD

Fig 3 (right): Brookfield DVE-LV Viscometer

Figure 4 (below): Disc spindle



Vertical Lift System

- Designed a dual scissor jack lifting mechanism connected directly to the electrospinner frame
- Lift is powered by a dual shaft DC motor coupled to lead screws and a potentiometer (controls speed)
- Roller wheels and aluminum t-slotted framing enable smooth and robust actuation

Viscometer Tests

- Viscosity is calculated from the torque required to rotate a spindle (Figs 3-4) in a solution at fixed speed
- Disk spindles give apparent viscosity, while cylindrical spindles give more accurate viscosity & shear values
- ω = spindle ang. velocity (rad/sec), R_c = beaker radius (cm), R_b = spindle radius (cm), x = radius of calculation (cm), M = torque (dyne-cm), L = effective spindle length (cm)

Shear rate (sec^{-1}): Shear stress ($\frac{dynes}{cm^2}$): Viscosity (Poise):

$$\gamma = \frac{2\omega R_c^2 R_b^2}{x^2(R_c^2 - R_b^2)} \quad \tau = \frac{M}{2\pi R_b^2 L} \quad \eta = \frac{\tau}{\gamma}$$



Fig 5: Upgraded solution delivery syringes (McMaster)




Fig 6: Channel adjustment model



Solution: PVDF, beaker diam. = 2", disc spindle LV-2

PVDF Weight %	RPM	Torque %	Viscosity (cP)
10	30	1.9	19.00
10	50	3.5	21.00
10	60	4.5	22.50
20	30	26.5	265.0
20	50	47.2	283.2
20	60	58.0	290.0

Results & Future Work

- Acquired funding and ordered components for the scissor lift design
- Procured new solution delivery syringe (Fig 5)
- Redesigned spinneret channel to include alignment adjustment set screws (Fig 6)
- Build the scissor lift and integrate it with the electrospinner
- Complete viscometer calibration and redo measurements of polymer solutions
- Incorporate new delivery syringe and channel adjusters

This manuscript has been authored by Fermi Research Alliance, LLC under Contract No. DE-AC02-07CH11359 with the U.S. Department of Energy, Office of Science, Office of High Energy Physics.
This work was supported in part by the U.S. Department of Energy, Office of Science, Office of Workforce Development for Teachers and Scientists (WSTS) under the Science Undergraduate Laboratory Internships Program (SULI).