




Mu2e Tracker Straws

Chiho Wang

Straws L3 Manager

10/22/2014



Outline

- Straw requirements
- Straw property measurements
- Straw assembly components and procedures
- Cost

Straw requirements


- 5 mm ID metalized Mylar tube cathode.
- Minimal material
 - 6 μm Mylar + 3 μm adhesive + 6 μm Mylar double helical wrap
 - Inner wall coating: 500 \AA Al + 200 \AA Au
 - Outer wall coating: 500 \AA Al
- Operates in vacuum
 - Sustains > 1 atm pressure difference
 - Leak rate < 7 ccm / detector volume
- Stability
 - Straw straightness: max. transverse deviation/sagging < 300 μm for HV stability.
 - Longitudinal tension is applied to keep straw straight.
 - Initial tension need to be higher to counter for material relaxation over time (creep).
 - Sustains radiation over the life time of operation.

Straw properties

- Destructive tests:
 - Pressure: sustained 60 psi (10 min.). Destroyed ≥ 70 psi
 - Stretch: sustained 1.6 kg (2 yrs.). Destroyed ≥ 2.9 kg
- Mechanical:
 - Linear density: 0.34 ± 0.01 g/100cm
 - Wall thickness: 15 μm
 - Derived from linear density, assuming Mylar density 1.39 g/cm^3 , and same for the polyester based adhesive.
 - Spring constant: 0.891 cm/kg/100cm
- Electrical:
 - Cathode resistance: 120 Ω /100cm

Straw leak measurement 1


Duke C. Wang, S. Oh



Pressurize to ~ 1 atm and measure the pressure drop

Straw leak measurement 1


Duke C. Wang, S. Oh



- N₂ leak rate = 0.095 ± 0.026 psi/day
- CO₂ leak rate = 0.36 ± 0.06 psi/day
- Assume Ar leak rate similar to N₂, and straw volume ~ 0.3 m³, estimated Ar(80%)CO(20%) leak rate ~ 2 ccm / tracker volume

Straw leak measurement 2




York J. Popp, K. Lynch




- Tested $10 \times 129\text{cm}$ straws
- Early results $\lesssim 4 \text{ ccm}$ / tracker volume
- Improvements underway

Straw leak measurement 3

Rice D. Rivera, M. Corcoran



EE891 and EE892 CO2 sensor




- Tested 119 × 129cm straws
- Averaged leak rate = $5 \cdot 10^{-5}$ ccm/straw
- ~1 ccm for tracker volume

Straw straightness measurement

To determine the required straw tension

Duke C. Wang, S. Oh

Horizontal Traveling Microscope (X)




Straw

Vertical Traveling Microscope (Y)

Straw deviation measured for different straw length and tension.

Straw straightness measurement


Duke C. Wang, S. Oh



Straw creep measurement


Duke C. Wang, S. Oh

- Glue straws on a support frame (120 cm) with tensions: 300 g, 400 g, 500 g, 600 g.
- Measure straw tension by resonant frequency as a function of time.



Straw creep measurement


Duke C. Wang, S. Oh



Corrected for temperature & humidity variation

Straw creep measurement

Duke C. Wang, S. Oh



Radiation aging study

Fermilab *R. Wagner*

- A single straw detector, constructed with relevant components, was operated under Sr⁹⁰ irradiation.
- Expected total dose (including beam flash) of 0.9 C/cm was irradiated.
- Gain change monitored by Fe⁵⁵ peak amplitude.
- Observed no measurable degradation in gain or cathode resistivity

Fractional gain change over 3 irradiation periods.
 Measured at 3 locations from downstream.
 Referenced to un-irradiated point at 110 cm:

| Period | Charge (mC/cm) | Current (nA/cm) | Gain(x)/Gain(110 cm) | | |
|--------|----------------|-----------------|----------------------|-----------|-----------|
| | | | p1 (7cm) | p2 (37cm) | p3 (67cm) |
| 1 | 120 | 18 | 0.979 | 0.982 | 0.994 |
| 2 | 120 | 36 | 1.010 | 1.009 | 1.002 |
| 3 | 670 | 70 | 0.994 | 0.990 | 0.983 |
| Total | 910 | | ±2% | | |

- Gas flow rate = 2 vol/hr/m
- 20±0.02% CO₂ with balance Ar
- Airgas “primary standard” grade. Purity 99.99% with individual impurities below:

| Contaminant | Max. Concentration (ppm) |
|--------------------|--------------------------|
| Carbon Monoxide | 10 |
| Hydrogen | 10 |
| Nitrogen | 20 |
| Oxygen | 10 |
| Total Halogens | 10 |
| Total Hydrocarbons | 10 |
| Water | 10 |

Straw property summary


- Straw is robust and has no problem operating in vacuum. The burst pressure exceeds 60 psi and is well above 14 psi.
- Straw leak rate of $\lesssim 4$ ccm is within acceptable range of 7 ccm
- 700 g initial tension is needed to keep straw tension above 250 g after 6 years, which is needed to keep straw straightness within ± 0.3 mm for HV stability.
- Straw can sustain tension > 1.6 Kg and is well above the needed 700 g initial tension.
- No measurable degradation in gain after expected dosage of radiation.

Straw assembly

- After a straw is received from manufacturer it needs to go through QC procedures and get assembled into a straw sub-assembly. These procedures include:
 - Visual inspection
 - Continuity check (measure resistance)
 - Leak test
 - Cut to length
 - Assemble/Attach end pieces
- A QA/QC data base is being developed to track/link component data.



Straw cutter

Fermilab *R. Wagner*




Straw sub-assembly

Fermilab *R. Wagner, A. Mukherjee*




Sense wire is aligned by external precision fixture (PAAS), not by straw components

- Simulation studies show it is *not* necessary to deaden straws near gas manifold



Old

- 4 piece construction
- Kapton sleeve trimmed to deaden ~1 cm past manifold
- Mostly brass





New

- 3piece construction
- Mostly plastic
- Lighter & cheaper

Assembly jigs for assembling end pieces


Fermilab *R. Wagner*



Assembly jigs for gluing end pieces to straw



- Orient end pieces on two ends
- Glue to straw with conductive glue

Fermilab *R. Wagner*



A 2x4 prototype

Fermilab
R. Wagner, V. Rusu



QA/QC data base

- Use a RDBMS
FNAL supports PostgreSQL
- Barcode users & most parts
Straws via their storage tubes
- QC and construction steps recorded in DB
- Leverage NOvA software experience and barcode readers, but write new software interface

Rice M. Corcoran
York K. Lynch



GUI

These would come from a barcode ...

The screenshot shows a software window titled 'Measure Straw Resistance'. It has a menu bar with 'File', 'Lookup', and 'Straw QC'. The main area contains several input fields: 'Straw' with the value 'Straw 2514', 'Operator' with 'Operator 2', 'Environment' with 'Temperature [C]' at 37 and 'RH [%]' at 32, and 'Resistance [Ohms]' at 132. At the bottom are three buttons: 'Measure', 'Confirm', and 'Cancel'.

... while these would be measured live from hardware

Program logic ensures consistency before insertion/update

Straw assembly summary



- Straw end components had been manufactured and tested.
- Assembly tools and test equipment for straw assembly procedures were developed and functional.
- Using the material and tools developed, a 2×4 prototype is completed with front end electronics attached.
- A QA/QC data base is being implemented for production process.

Cost Table

| | Base Cost (AY K\$) | | | Estimate Uncertainty (on remaining costs) | % Contingency on ETC | Total Cost |
|-----------------------------|--------------------|-------|-------|--|----------------------------|---------------|
| | M&S | Labor | Total | | | |
| 475.06 Tracker | | | | | | |
| 475.06.02 Straws | | | | | | |
| 475.06.02.01 Straw Tubes | 1,170 | 67 | 1,237 | 425 | 38% | 1,662 |
| 475.06.02.03 Wire Stringing | 29 | 2 | 31 | 4 | 21% | 35 |
| Grand Total | 1,200 | 69 | 1,268 | 429 | 38% | 1,697 |


Cost Breakdown

Base Cost (AY K\$)



Quality of Estimate


Base Cost (AY K\$)



- L1 Actual / M1 Existing P.O.
- L3 / M3 Advanced
- L4 / M4 Preliminary
- L5 / M5 Conceptual


Labor & Material by FY

Base Cost (AY K\$)



Labor Resources by FY

FTEs by Discipline



Conclusion

- Detecting elements of the tracker: straws and straw-end components are manufactured and tested.
- Straw tube is robust and leak tight, and exceeds the requirements to operate in vacuum for 6 years.
- Assembly procedures and assembly tooling are realized and exercised.
- A QA/QC data base is being implemented.
- Ready for CD-2